

Product datasheet for **SC213763**

EVI5L (NM_001159944) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EVI5L (NM_001159944) Human 3' UTR Clone
Symbol:	EVI5L
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001159944
Insert Size:	1304 bp
Insert Sequence:	>SC213763 3'UTR clone of NM_001159944 The sequence shown below is from the reference sequence of NM_001159944. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCGCCCTACAGCCAGGGTCTGGACAACTAGGGCCATGCCAGCGCGCCCGGAGTCAGGAGGCCGAGCC
GCGGGGGCGCCCGGGCAGTCCGCGTTCTGCTCCCCACCTGCCCACTTGACAACTACGCGCCCTCTG
TGGCTCGGCCACCCCTAAAGCGAGGCCCGGCGAGGCAGCGCAGAGGGTAGGGTCCGACCTGGGCTCCTC
AGGGCCCGGGGAGGCTCTCTATCCCCAGCAGTGTTCACCATCTTGGTCTGTACCCCTCCGGGCCCT
CTGGCGTTCCAGGGGTGCCTGGAGGGGCTGACTGCTCTCTTAACAGGAGGGCAGAGGGCAGGGGACAGA
CGACCCAGAGGTCCCAGCACTGAATGAGCAGGCAGCTCCACCTCCTGGCAGGCTTCCTTCTGGGACAG
GGGCGACATTGCTGGGAAGTGCTCAGGAGGTAGCCGAGGCCTGAGGAAGGAGAGCGCCAGCTCTGGGCT
GGACATCAGCACCCACAACTCCTCGGGATGAAGTGACCCCTTACTAGCCCTGGCCATCTCTAGGG
GAGTCAGGCCGCTGGGGACAGATGGCCAGGCCGGCCTCTCCTGCCTGGCGCAGGCACCGTGGCCCTGC
AGCGGAAACCAAAGTCCCCCTACTATGTGCGGGGCGCCTGGGGACTGAGTGGCTAACGGGAAGCCTCCC
TGCTTCTCTGGGGCCAGAGCAGCTTCAGGAAGTGAAGAGACCCCTCTCCCCAGTGCCCCACATCTTTC
TCTGGGAGACACTCGCGCCCCCTATCCTGGCCAGTGATGGAGGGGTATCTCTACAGGGGTCTCGGTCT
CCATCTACTTCCCTTTCTTCATGTGTTTCCAGCCCCACCCTCAGCCAGAGCCAGGCCCCAGGCAGGAG
CTTCCCAGCGAGCGGCCTCCCCTCACTTCTCCTGGTGGTCCCCTTGTCTCTGCTGAATCAGAGCTGAG
AACGGTGCCAAAATCCAACGATGCCCCCAGGCCCTTCTCCTTCCCCCGGCCCGGGGCTCATG
ACCCTGGGGGCTGCCCCCTCGGTGCTCCCGGAGCCTCCAGTAGAGTAGCGTCACAAGCAATCTCCCTG
GCGCTTCTGGGTGGGGACCCCTGCTCCGTCCCGCTTCTAGCTGCCCACTTTTCAGTGTTACGAAGC
CTGGGGACCGGGCAGGCACCCACGGGGCTCTCCACACGCCCCCTACACTGCCCGCCACCATTTTGAC
ACTGCCTGTTACATGTCGCCCAGGCGGGAAAAATGAAAAATAAGTGTATTTACACAGTCA
ACGGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



View online »

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001159944.3</u>
Summary:	Functions as a GTPase-activating protein (GAP) with a broad specificity.[UniProtKB/Swiss-Prot Function]
Locus ID:	115704
MW:	46.1