

## Product datasheet for **SC213765**

### **EVI5L (NM\_145245) Human 3' UTR Clone**

#### **Product data:**

|                           |                                                                                                                                                                                                                                                    |
|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type:             | 3' UTR Clones                                                                                                                                                                                                                                      |
| Product Name:             | EVI5L (NM_145245) Human 3' UTR Clone                                                                                                                                                                                                               |
| Symbol:                   | EVI5L                                                                                                                                                                                                                                              |
| Mammalian Cell Selection: | Neomycin                                                                                                                                                                                                                                           |
| Vector:                   | pMirTarget (PS100062)                                                                                                                                                                                                                              |
| ACCN:                     | NM_145245                                                                                                                                                                                                                                          |
| Insert Size:              | 1304 bp                                                                                                                                                                                                                                            |
| Insert Sequence:          | <p>&gt;SC213765 3'UTR clone of NM_145245</p> <p>The sequence shown below is from the reference sequence of NM_145245. The complete sequence of this clone may contain minor differences, such as SNPs.</p> <p>Blue=Stop Codon Red=Cloning site</p> |

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GCGCCCTACAGCCAGGGTCTGGACAACTAGGGCCATGCCAGCGCGCCCGGAGTCAGGAGGCCGAGCC
GCGGGGGCGCCCGGGCAGTCCGCGTTCTGCTCCCCACCTGCCCACTTGACAACTACGCGCCCTCTG
TGGCTCGGCCACCCCTAAAGCGAGGCCCGGCGAGGCAGCGCAGAGGGTAGGGTCCGACCTGGGCTCCTC
AGGGCCCCGGGGCAGGCTCTCTATCCCCAGCAGTGTTCACCATCTTGGTCTGTACCCCTCCGGGCCCT
CTGGCGTTCCAGGGGTGCCTGGAGGGGCTGACTGCTCTCTTAACAGGAGGGCAGAGGGCAGGGGACAGA
CGACCCAGAGGTCCCAGCACTGAATGAGCAGGCAGCTCCACCTCCTGGCAGGCTTCCTTCTGGGACAG
GGGCGACATTGCTGGGAAGTGCTCAGGAGGTAGCCGAGGCCTGAGGAAGGAGAGCGCCAGCTCTGGGCT
GGACATCAGCACCCACAACACTCCTCGGGATGAAGTGACCCCTTACTAGCCCTGGCCATCTCTAGGG
GAGTCAGGCCGCTGGGGACAGATGGCCAGGCCGGCCTCTCCTGCCTGGCGCAGGCACCGTGGCCCCCTGC
AGCGGAAACCAAGTCCCCCTACTATGTGCGGGGCGCCTGGGGACTGAGTGGCTAACGGGAAGCCTCCC
TGCTTCTCTGGGGCCAGAGCAGCTTCCAGGAAGTGAAGAGACCCCTCTCCCCAGTGCCCCACATCTTTC
TCTGGGAGACACTCGCGCCCCCTATCCTGGCCAGTGATGGAGGGGTATCTCTACAGGGGTCTCGGTCT
CCATCTACTTCCCTTTCTTCATGTGTTTCCAGCCCCACCCTCAGCCAGAGCCAGGCCCCAGGCAGGAG
CTTCCCAGCGAGCGGCCTCCCCTCACTTCTCCTGGTGGTCCCCTTGTCTCTGCTGAATCAGAGCTGAG
AACGGTGCCAAAAATCCAACGATGCCCCCAGGCCCTTCTCCTTCCCCCGGCCCCCGGGCCTCATG
ACCCTGGGGGCTGCCCCCTCGGTGCTCCCGGAGCCTCCAGTAGAGTAGCGTCACAAGCAATCTCCCTG
GCGCTTCTGGGTGGGGACCCCTGCTCCGTCCCGCTTCTAGCTGCCCACTTTTCAGTGTTACGAAGC
CTGGGGACCGGGCAGGCACCCACGGGGCTCTCCACACGCCCCCTACACTGCCCGCCACCATTTTGAC
ACTGCCTGTTACATGTCGCCCAGGCGGGAAAAATGAAAAATAAGTGTATTTACACAGTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

|                           |                                                                                                                                                                                                                                                                                            |
|---------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Restriction Sites:</b> | Sgfl-Mlul                                                                                                                                                                                                                                                                                  |
| <b>OTI Disclaimer:</b>    | Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs). |
| <b>Components:</b>        | The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.                                                                                             |
| <b>RefSeq:</b>            | <a href="#">NM_145245.5</a>                                                                                                                                                                                                                                                                |
| <b>Summary:</b>           | Functions as a GTPase-activating protein (GAP) with a broad specificity.[UniProtKB/Swiss-Prot Function]                                                                                                                                                                                    |
| <b>Locus ID:</b>          | 115704                                                                                                                                                                                                                                                                                     |
| <b>MW:</b>                | 46.1                                                                                                                                                                                                                                                                                       |