

Product datasheet for **SC213772**

p90 Autoantigen (KIAA1524) (NM_020890) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	p90 Autoantigen (KIAA1524) (NM_020890) Human 3' UTR Clone
Symbol:	p90 Autoantigen
Synonyms:	KIAA1524; p90
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020890
Insert Size:	1315 bp



[View online »](#)

Insert Sequence:	>SC213772 3'UTR clone of NM_020890 The sequence shown below is from the reference sequence of NM_020890. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site
	<pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AATCCAGAACTGTGAATCTCAGTATAAGACATTATGGCATTGGAATTTGTAATCTCATGATATTT TTGATGTATTTATCTATTGGAGGGGGGTGGGTAGGGGAGTTAATTTGTGACTTCGTAACAATAAGAAG TTATTATCTAATTTAGTAAAGACCCTGATCTGTTGCATGTTTTTATTTGATAGTTTGAATAGAAATTT AATTTTCTAAGTTTTACTTTTTGTTTCTGGCTTTTATGGCTTAAGGTTTTCTTTGGTTCTTACATTAGA AAATCATTTTTAACCTCCATTATCATTTTTCTAAGGTTCTTCTTTTCTTAGTTGCTTCTATTCTG TTTTGCCTGTCTATTTATTTCTCATTTGTGATTATTTAGATCTTAAGACCAAACCTTCTTGGTATAACA GTCCTAAAGATTACAAAATAAAAATATAGAGAGAGTAAAGTAAAGTAAAGTAAAGTAAAGAGAGAAGGT ATGATTTTACTGCTAGAGAAGTTGTCTCTGAAGAAGCACAAAAGAAAATATTAGTGAATTTAAATAA TTTTTATACTGCTGTAGCATAATTTCTAAATTTGAAAAAATGCAATGGTAATAAATGTATAAAAT AGAAAACGTGCTATTGTGTAAACTATTACATTTAAATGATTACATTTAACACAATAGCTGTCTCATAAA AAATCTAAGAACTGTAGAATTATTTGTAAGGTATATTTAGTGTTTTTCTCTTTTTTCAAATGTTGC ATAGTGGTCCAGATCATCTATAATTAATTATAAGTTTCTTATGATTTACATAAGGATAATTTGGATTCT ACTCTTGCTGCATATTCAGAAAAATATTATTTACTGTTAATTTCTGATATGAATTGATTTAATTTGTTG TCCTAATAGTTTTTAATTGGTTCATTTAATTTAAAGCATAGGTTAAATATCTTTATTTTTTTTAAAT ACAAGGGTCCTTCTGACTTTTCATTTGAATCATATCCTCTAACTACTCTGTAAGAAACGTTTTTACAT CAATATTTTGATTGATTTTATTGCTATACCAATACAGTCATTAATTTTACTTGGTGTTTTCTTCACCAT ATGTCGTGTTTGATTTAAGTGAATTTAGATTTGTACAGTTCTGGAATTTTTTAAAAAATATTATGCAGA GTAAACATTTTAAATTAATTTGTTCTATAGAGGAGATACATTCTGTTTCTCTTTTGAATTGCATG CATTTTCAACCAATATGAAAAAGTTTAAAGTGGATAGATAAGACTGAAGAATAAAAGGCTATAAAAT AAAA ACGCGTAAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_020890.3
Summary:	Oncoprotein that inhibits PP2A and stabilizes MYC in human malignancies. Promotes anchorage-independent cell growth and tumor formation.[UniProtKB/Swiss-Prot Function]
Locus ID:	57650
MW:	52.2