

Product datasheet for **SC213780**

UBE2J2 (NM_058167) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: UBE2J2 (NM_058167) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: UBE2J2
Synonyms: NCUBE-2; NCUBE2; PRO2121
ACCN: NM_058167
Insert Size: 1321 bp
Insert Sequence: >SC213780 3'UTR clone of NM_058167
 The sequence shown below is from the reference sequence of NM_058167. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TACGTGCTGAGGAGCATCGCGCAGGAGTGAAGGCCAGGCGCCGAGACCCAAGGCCCACTGAGGGCACC
GGCACCAGAGCGTGACCTCGGCAGGCTGGACACACTGCCCAGCAGAGCCAGACCACCAGGCTCCTAG
GTTTAGCTTTTAAAAACCTGAAAGGGGAAGCAAAAACCAAAATGTGTACTGGGCTTGGAGGAGACTG
GAGCCTCAGCCCTGTCCTGGCCACGGGCGCTGGGGCTGGTGTGGGTGGGCTTGTGTGCTGGATTGT
AGCTTATCTTCCGTGTTGTCTTTGGACCTGTTTTAGTAAACCCGTTTTTCATTTTATTAGATGTGGTCA
CTTAGAAATGCAAACCTGCTGCCAGCCGCGGCTGCTCCTGCGTCTTGGAGCTCCTGGCGGTTTCTC
GGAGCTCCCGGCTCCTCAGCGGGTGGGAACCTCGGGGCCAGGGGTGGAGCTGGCGTCCGCGGGTGTG
GTCTGGCCTGGCCGTGTGGTGATGAGGCTTAGCGGGCCAGTGACGGCCGTGGCTCAGGATCCATAAGT
CGGGGTTTGGTCTCAGCATTTACAAATGTGTTTACAGTCAGAATGAAACACATTCCTTCTAGAAAAGTGC
TTGGGGTTTTTGTGCCCTGGAAGCCAGGAGCCTGCTACTCCAACCACAAGTCGCCCTTGACTCGGG
CGGCCGCGAGCGGGCGGGGCTGCCGGTGCCCTCCGCGAGCCGGGCTCCTGGGCGCCCTCGGTGCT
GCAGGCTGGGGGCTTGGGTACCTGCAGAGCCTTTTCTGAATTCCTTATGTCCGGTGGGCCAGAAG
CCCGTCTCCTATGCTGGTGAAGGCGGAGGACCGGAGTCCCTGCAGAAGGCCCGTGCCTCGGGGGC
CTCCCTCATATCCCGTGCCTTGCCTGCGCTTACAGTAGGTAATGGCTCCGGCCCGGGTGTTCGCT
GTCCACGGAACATGGCAGAGGGGACCCCGGCCGAAAGACGCCAGAGCCAGCAGGGGCTGTTTCGGG
CCGCGTGGCTCCCGGGTCTCGGCCGTCTCCCTCTTCTGCGTGTGTTCCGTGACTTCGCTGGGTGGG
ATGTACCGCAGGTGCATCGCGTCGAGGTGGGGCACGGCCCGCCGCAAGAAACCCACCCTGTCGGAGGC
GGGCGTGAGACAAGCCCAGCCCAGCGCCTCATCTTTCTCGTTTTTTGATCAGTTTATTCAGAATTG
CTCTATAATTTACCAATTGTATGATTTAACCTATTCTTGTGAAAAAAAAGGTCTTTCATTATATCTT
TATTTCTGAA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_058167.3
Summary:	The modification of proteins with ubiquitin is an important cellular mechanism for targeting abnormal or short-lived proteins for degradation. Ubiquitination involves at least three classes of enzymes: ubiquitin-activating enzymes, or E1s, ubiquitin-conjugating enzymes, or E2s, and ubiquitin-protein ligases, or E3s. This gene encodes a member of the E2 ubiquitin-conjugating enzyme family. This enzyme is located in the membrane of the endoplasmic reticulum. Multiple alternatively spliced transcript variants have been found for this gene, but the full-length nature of some variants has not been defined. [provided by RefSeq, Jul 2008]
Locus ID:	118424
MW:	47