

Product datasheet for **SC213816**

IPO11 (NM_016338) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	IPO11 (NM_016338) Human 3' UTR Clone
Symbol:	IPO11
Synonyms:	RanBP11
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_016338
Insert Size:	1294 bp



[View online »](#)

Insert Sequence:	>SC213816 3'UTR clone of NM_016338 The sequence shown below is from the reference sequence of NM_016338. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CAGCTACAGGAGTTTTTGCAGGATTCTAAGAGCACATGACATGTGGCTGCCTCCCTTTTCAGAAACAA GCTGAGTAACCCAGCCTGCCGTTTGTATGTGAGAGCCTGCTGAGATGAAGAAATCACTTCATGAAAATA AGCAAAGACCACACATTTTTTACTACAAAATGTAAAGGATAAATGTAAATCCTGCATACTAAAATCAC AAACCTATTCTCAAAAGAATTTAATTTTATATTTATGAGGGGGCCCTTCACTAAAAAGTACATGTAAA AGTACATTTGATGACAATAGCTGCTTAGTTTCTGTAAAGAGAAGAACTTTATCTTTAATTATGTGC TCTTAATATTTGAAGATGAGAGTTAATACCTGAGATGTTTTCTGCAACCAAAATTCATTAAATTTGGC TGCCTTATCCTTTTTTTAAGCTAATGAACTACAGGTTTGAAAAATGACAAAGCTGTTTCAGATGATGCT ATTAAAGAAATGTGTGTAAGCAAAAATATATAAATAGTGACAAATACACATTACCAAGCTTATCTT GCAAGGGAGTTATTTTCATCTAACATAGAAAGTGTTTTATCAGACAAATGCTTTTATTTTCATTCTA ATAATTTGATACAGAAATTAGTAAAGGCATTTTTTTCTTTTTTTTCCAGTAAATACATTGGGTCTATA AATGTGCATTTGTAAGGGCCACAAAAGTGAACGTGTGGTACTGTAGTACCACGTGGGAGACCTCTGGTT ATGGTTTAGTCCTAGTTCCTTTGTTACTCCTGTGAGCACCGAGAAGAACTGGGCGACTCCAGTCCCAC CTGTGCTGTGACAGTCCCACGTGGCTATGACAGACTGTTAGTACTTACCCTTCTCAGGTTCTCAGTG CAGGGGTGCATCAGGGCCTCAATAATAGGGGTATACCTGGGAGGATCCAGCAGTAATCCCAGGGTACT AGGATTACTAGTACTCTGATGGAAGTACTTCTTCTTCTTATTCCTCGAACATGCAGTACATAAAAAGG GGAAAAGGAGAAAAAAAAGCCTTACTTTGTTTTACTTGCCATTTATTGTAAGGAACTTTAAAGCATT TTTTAGGAAATACTCAAAAGCAAGGTTGGAAAATGTTTTATCTTTCTATAGAAAGTTGGGTACAGTATG TAAGTGCAGGAAACCACTGCCCTTTGTAAGCTGTGGAACCAAACTGTATGGGATATTTGATGTTT TCAGAAAGAGGAAGAAAATATGGTCCAAATTAATTTTCCAAAGATACCTCA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_016338.5
Summary:	Importins, including IPO11, are a members of the karyopherin/importin-beta family of transport receptors (see KPNB1; 602738) that mediate nucleocytoplasmic transport of protein and RNA cargoes (Plafker and Macara, 2000 [PubMed 11032817]).[supplied by OMIM, Sep 2008]
Locus ID:	51194
MW:	49.7