

Product datasheet for **SC213845**

LRP3 (NM_002333) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LRP3 (NM_002333) Human 3' UTR Clone
Symbol:	LRP3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_002333
Insert Size:	1319 bp



[View online »](#)

Insert Sequence:	>SC213845 3'UTR clone of NM_002333 The sequence shown below is from the reference sequence of NM_002333. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGCGATCGCC AGCGATGATGAGGCCCTGTTGGTCTGTGACCGCTGGGCTCGCTGGTGACCGCCACAGCCCCGCTTTGT AACCAGGGAATACACAGTCATTTCTACCTGCCTCTGCGTCCTTTCTTATGGAGAGGCCCTCCGGGGAC CCCAGCGGAGGGGCTGGCCCTAAGCCAGCTGGCTGCACTGGTGGGCGGGAGCTGTGGGACTGAACGGC GGGGGGGAGAAGAGTGGAGTGGTGAAGCCGCTGCAAGGGTCCCATTGTACACAAGCACCCTGGGGTCTC ACTTCTCTCCCCCACTCCATTCTGGGAACCCATTTCCAGAGAAGCAGGGGACCAGAGCCTTTTTGCTT CATCTGCCCTGCAGTGGCAACAGCTGCCCTAGAGCTGGGAGTTGTCATGAGGATGGGGCCGGCCATG GGGACGCTGGGTCTCATCCGTGGTGACTATTTTGCTCACGCCTCACCTCTTCTGTCCAGCAAGGCCT CTGGAGGACCTGGGTGGGTGGTCTCTGCCAAACCCTCATGCCCTGGCCAGGCAGCGCCCTTTGCC TCCTTGAAGGTCCAATCTCCTGGGCACCGTGGAAGGGCTGGGGTGTGTGCTGTGTACATGAAGCCCC TGGAGTTCCCTGCAGGCCTCAAAGCCCCCAGGACATGCTGCTCTGTGCTTTCTGGGGACAGAGAGCCG AGGCACCAACAGCAGTGGCCACGTCTGTTCCATCCTATCTCAGGGCCTCGGTTTCCCCACCTGTAAGCC AGGGTGAGACTGTGGGCTCTGAAGTCTGACCTTCACTCTGTGGTGATTACGGGAAGAAAGACAATCC CAGCCTGGAGCCCTGAGCTCTCCTGGGTGGCCGAGCTCAGACTGGCATCGAGGGGCCTGGGGTAGGGG CGGACTGAGCTGTACCCCATCTCGACATCCCTTTAGGGCAGGGGGACTCTGGGTGCCATGGTAGGGA GCGCCTGTGTCAAGCAGCAAAGCCCCCTAATGTCAAGAGCCTCTCCAGGGCCTGCCACGGCATCTTCC TGCCAGCAACATGTGTGCACCCACCGAGGTACGCCCCAGCCACTCCCATCTGTGCCCATCAGGGACTC CCATAGCACGAGCGAAGCAGCCAGCCCTGTTTATTTATAGGCCTTTTCAGGAAGAGCTAGCAGGGCAGT GCTAAGACAGGAAACCCAGTCCACATTTTAGGGCTTCTTAAACAGGCTTCTGAGAGTCGTATCTTTT TCTTTTTTTTCCAGAAAAAACAAAACTTTTTTGCCAAACACCTCTCAATAAACACATGTA AACAGAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_002333.4
Summary:	Probable receptor, which may be involved in the internalization of lipophilic molecules and/or signal transduction. Its precise role is however unclear, since it does not bind to very low density lipoprotein (VLDL) or to LRPAP1 in vitro.[UniProtKB/Swiss-Prot Function]
Locus ID:	4037
MW:	47.1