

Product datasheet for **SC213866**

COL24A1 (NM_152890) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	COL24A1 (NM_152890) Human 3' UTR Clone
Symbol:	COL24A1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152890
Insert Size:	1343 bp



[View online »](#)

Insert Sequence:	>SC213866 3'UTR clone of NM_152890 The sequence shown below is from the reference sequence of NM_152890. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC ATTGACAGCAGTTCTGTATGCTTTCTGTAAGTCTCTGAATTAGTTCCGAATTCAGGCTGTTGGCCAGG TAATTGCTGCAGAGGGAGAAATAAGACAGACAGATACAGTCATTATGAAATGCATGTAATAAAGCATTG GCTAAATCTTAAAGAATCTCAGGAAGAACAGACTTCTCCTAAGAAGGAGAAAAGGCATTTTAAAGGA CTATGATTGATAAAGTATTTAATCTTTTAAAAATTATATTCATCTCAGCTTTCTTAGAGAATCCCTA GAACTAAAAATTTATAAATATGGAATCTTCAGGGTATCTTATTTTTGACTGAGTGCAGTACCCA TTAGACAGCTGGAGATGCAGAGCACTATGGAGCAATACTGGCTAATGCTTCCAGATGTGCACTGCTTCT GTCTAAAAATTACAAGCCACAGTCTAATATGTCTTATTTTCCAAACACTAAGCTGTATTCAGGTCCCC GATGGGCATATACATCTTAGCCGGTGATACACTACCTCTTACGTGTTGCCTCTTTGTGTTGCTTGGTGC TCTTTTCAAAAACAAGGTGCTTATGGCTTTCATAGACTATTTCTTTTTCATCTTTGTCTTTTAAAA GTGTATGTACTGGTTACATCAAGATATGTTTTGGTTGTTAGTACTATTTTAATTTGTTTGGTCACACA CTTAATAACACGTGAACTATTTATGTGAAGTCCTTGTGTTTTATTTTAAAAATCTCTTTGTGATTGGA ATCAAAGCCAGCAGATTGTAACCTGTGCTTGTACGCAAAAGAAATTAGATTTCTTTGTTTTGTTTTATT TTTTAAATTGTTGTAATAATTATTATAGGCCAGCTACATCTAGTAGTAGTTTGGGGTACAGATTGGGG GTTGTGCCACTGTTTTTAAAGTTCATGATCATCTGGAATGATACTTAGTGATATATATTTTGTAAA GTTTTAATTCAGCAAATTTTTTAAAGTTCATGCTGTTTTTAAATTATAAAACCTTTATATTTCTGCTTTG TAGAAATTATATGTTTTGTAGTATTCATTGATTTTCTTCACTGTACTTAAATTTAGTGTTAGTACTTT AAAATTTTTTAATTTACCAGTCTTTAAAGCAACATCCAGAAAAAAAAGTCTTTTCCCATTTAAAAATAG GCTCAGCCAGTTCAATGTCTCCTTGTTATCAGAGAAATATTAGTTCAATACTGAAAGAAAAATATTATA CCTCTTGGTATCTAGAAAAGCTTGTTCATCCATTATAATATATCTTTAGCCACAGCAAACCACTTA ACCTATCTATAATAAAAAATGTGCTTTAAATAA ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152890.7
Summary:	This gene is a member of the collagen gene family and is thought to regulate type I collagen fibrillogenesis during fetal development. [provided by RefSeq, Mar 2017]
Locus ID:	255631
MW:	51.8