

Product datasheet for **SC213881**

NCOA5 (NM_020967) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	NCOA5 (NM_020967) Human 3' UTR Clone
Symbol:	NCOA5
Synonyms:	bA465L10.6; CIA
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020967
Insert Size:	1339 bp



[View online »](#)

Insert Sequence:	>SC213881 3'UTR clone of NM_020967 The sequence shown below is from the reference sequence of NM_020967. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CCCATGGGATCTTACCAGAGGCATTACGAAGCTAAATCTTTCAACTCTCCCAGTCCCCTCATCCCCT GGCTCTCTCCCACTTACTTGTCTAAATAGAGCTGTTTGAGGATGTTCTCTGCGCTCCCAGGCCGCA TCGAGTGTCAATTTCTACCACCTGCTCTCTTCTGCCAAGGCTGTGTGTTATTCTTACAAA GTTTATACTGCATTTGGGGCTGTATCTTTTTTGTGTTTTGTGAGAAAATAAATATCTCCGGGG CAGTACAGGTGTCTGGGCTGTATTTGATGGGTTTCTCCGGTCCCTGTTTCTACTGGATTGGAGCCA GGCCAGCTAGCCAAGTTTGAATGGCATTGTGTCATGTCAGTAGCCACCACCTTTGTTTCATTGTGAACC TACCAAGGCTTTCCAGCTTCATACACATTGACCAGAGCTCAAGCTCCTGCCTGCAACTCCTGCCTAGAG TTGAAGAAAAGCAAAGTGGCCTTGGCAGGCACAGTGTATCATACCTCACCCCATATGTTTGGGTCT GCTTGAGGATTATAAATCAGCCACTCTGGATTGTTGAGGAATGGCCATGGCAGCCACAGAAAAAGAA TTTTTCTCTGAGCCAAGGTTGTTTTTGTGTTTTTCTTTTTCTTTTTTGTGTTTCATTTTCATTGAA GATCTCCAATGGACTGAACAGCTCCAGTCAGCAGCAGTTACCACAACTGTGAATCTGGGCCCCACCAC TCTTCCCTGTTAACCAGTTCTGTGAGCATCCCCCTCTCCAGCAGCACTTCCATGAAGTTGGTTCTGAGA CTCTGGCCGTGAACACCCGTTTCTTCACTGATTGTTTGGGCTTTGGCTCAAAACCCAGGCTCTTG TTTTGTCTAGACTCTTATTCTGTTTCTGAGCAGCAGGAGGTAGGACCACTTTGATGTCAGACTTCT GGTAGCTGGACATGTTCTCGAGATGGGTGGCTGTTGCGCACTTTGTACCAGAGTGAAATGTTAGAAG GAGGGTTTCTGGCTGTGGTTCTAAATGGAGCCCCAGGAAGCTGCCCTCTCCCAGGTTTGTGCTCAGT AGAGCCCTGTGGATCACAGTCTTGAGGTCCTCTAGCAGGGGTGAGGGAGAGCAGCGACTTCAGCTGAGT CCCTGCCAGTGGTTAAGCAAACAATGGTTTCAAAATCAAGGTCCCAATGGCAGCATTTTATGTTCT GACCTGTTGTGTATATAGTGGTTTTTTTTTCTCTTGGAACTCTTGTTGTTAATAAAATGAAA TGATTACTTTTTAATTAATGAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_020967.3</u>
Summary:	This gene encodes a coregulator for the alpha and beta estrogen receptors and the orphan nuclear receptor NR1D2. The protein localizes to the nucleus, and is thought to have both coactivator and corepressor functions. Its interaction with nuclear receptors is independent of the AF2 domain on the receptors, which is known to regulate interaction with other coreceptors. Several alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jan 2017]
Locus ID:	57727

MW: 50.3