

Product datasheet for **SC213892**

Bestrophin 3 (BEST3) (NM_032735) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Bestrophin 3 (BEST3) (NM_032735) Human 3' UTR Clone
Symbol:	Bestrophin 3
Synonyms:	VMD2L3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032735
Insert Size:	1328 bp



[View online »](#)

Insert Sequence: >SC213892 3'UTR clone of NM_032735
 The sequence shown below is from the reference sequence of NM_032735. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACAAAGGAAACTGAGGAATCACCCAAATGAGTGCCACCAAGTTCTAGGACCTGGTCTAGCCTAGATTC
TTACCTCCCAAAAGCAGCATATTAGGCCCAGGACATACCTGAAGGCTGGTTCGACTTTTTAAAAACAT
GATCACACTGTATAAGCTACTTAGAACTTTGAAGGGCATTATGATCCTGACTTTTTGAGAGCTGCGAAA
AATCAAACGCATTACGTCATCCAACACAACATGGCTTTCATCAGAAGTTACACAGACCACAGTGATTT
AGCCAAAATAAGTCTTAACCAAAAAAAGGAATCTGAAAAACATACATGACTTATAAATTCATAGATTATT
AGTTGTGGAAGAGTCTCATAGATCATAGTCCAATCCCCTAGACCCAAAGGCCAGACACGCAGGGTTC
TCTATTTCTCTTTTTCTTTCTTTCTTTCTTTAGCCATAAAGAAAGGAGATGGATCACCTAATG
CCTGAAAGATCCCCTCCTCCTAACGTTTTTGTGGAGGCTTTCTAACCATCTCATTCTACTTTTATA
ACTATTATTGTTATTATCACTTCTGATTCTTAAGCTTTCATCAGCACCAGAACACTTAATTTCTCTCAA
TTTCCTTTTTTACTTTTTCTCTTTTTCTTTTATTCTGTGCATCCTTTTTTCCCTTCCCTTCTCCCTGG
TTGTGCTGCTTTTCTAGTAATCTGTCCCATTTTTCAGGTTATAATATTTCTCTTGCCAAATCTGTTGC
TTATGTCTTGCCATGAATTTCAAATAAAGACATTGGCAAAGCCTCCTACAGAATCATCAAAGGTGTC
AACAAATCCAGCCTGATACCCCTCCAGTCTCACTTCCAGTGGAAGGCATATAGGAGTGATGGTCCCCC
CGTGCCCTTGAGCCATTGACATTTCTTCTGCTTCCACTATTTGGTGTGGGACTGGAGAAAAACGA
AGCCAAGAGAATGAAATTTGCATAGCCAAATGCTATGGAACCAAGAATAAGATGACTCACATCTGT
TTAACTTGAACCTTATGCAAAAATCCTAATCGTGTGCAACACGGGACCAGTTACTGTGTCTATCGGG
TCATTTGGGAATGGGACGATACTGGAGCTGGGAGGCTAATGGGTCCAGAGCTGCATAATGGGGATCG
GATGTTGGTACTCTCCAGGGAGCTCTCAGGACAGCTGTTTGTGCTGGGCTCCCAAGTGCCTAGGCCA
TGCCTTACGACTGCTGCACTTACCATTTGCATATGATGGCACTTACTACATATGGCTTTAGAAATAA
AATAACCCTCAAACCTA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_032735.3](#)

Summary:

BEST3 belongs to the bestrophin family of anion channels, which includes BEST1 (MIM 607854), the gene mutant in vitelliform macular dystrophy (VMD; MIM 153700), and 2 other BEST1-like genes, BEST2 (MIM 607335) and BEST4 (MIM 607336). Bestrophins are transmembrane (TM) proteins that share a homology region containing a high content of aromatic residues, including an invariant arg-phe-pro (RFP) motif. The bestrophin genes share a conserved gene structure, with almost identical sizes of the 8 RFP-TM domain-encoding exons and highly conserved exon-intron boundaries. Each of the 4 bestrophin genes has a unique 3-prime end of variable length (Stohr et al., 2002 [PubMed 12032738]; Tsunenari et al., 2003 [PubMed 12907679]).[supplied by OMIM, Mar 2008]

Locus ID:

144453

MW:

50.2