

Product datasheet for **SC213917**

TSHZ1 (NM_005786) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TSHZ1 (NM_005786) Human 3' UTR Clone
Symbol:	TSHZ1
Synonyms:	CAA; NY-CO-33; SDCCAG33; TSH1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_005786
Insert Size:	1333 bp



[View online »](#)

Insert Sequence:	>SC213917 3'UTR clone of NM_005786 The sequence shown below is from the reference sequence of NM_005786. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ATCTATGTGACTGAGTTGGAGAAACAGTACGTCAGGTATGCAAGAGACCGCGGAACATTGCACTAAA CGTCGTCGAGCTGCACTAGGCCTGGCCTGAGCCTCTGAAATCAGTCTTTCCTTTGTTGCTGGCCCGCCT CTCTGGACCTTGGTTTTCTTACACATATTTGTATATTTATATGCTCTGTCCGATCTGTGCATGTTA TTTTTCTTTTCCGTGAGTCAAAGTCTGACCTTTATTTTCAACATCTGTTCTTGGTGTTAAGCTATCTT TTGTAGGAAATAGTGGGGCACACTACTCAGAGACATTATTTAGCAGTAAAGAAAGACACAAATAACAAT GATAAAAAGACATCCTAAAATGGTGAAGTTGCCATGACAATAAAGGTCATAGAACCTGGTAGTGCAAA TTTAACCCCTTTGAGGACTGTAATTGCATTTCTGTGCCTTTCACCTTGAAAAAAAAAAGAAGAAAAAA GGCTTAAAGGCATTTATTTCAGATCAAAAGCTGTGTCAGACACAAAACCTATTTAATAATTAATAATG AGCTTTTATATGCAGAACTGGTTTGTGAGGAAACATGCATGCAGCTGTGGGAAGATAGAAAGCTGTCTA GGACTCACGGGGTTTCGGAAGCCACACCTCTACATGTCTAAAAACAAAAACAAACAACAGATACCCCA CTCGCCACTATGCCAAACCTCTGGATGCTGAAAAATGCCACTTTCGGCAAAAAAAGTATGCAAG CTATTATTTAAAAATGTGTAAAAATGCTATTTCTTTTCTGTAAAAATATTGCAGTTTATAGTTATTT ACAAATGTAAGCTTTGGTACAAGCTGACCTTCTATAGCGTGTGCATTGGTAAATGAAAAAATGGG GCAACGTTGGAGTCCTTTCCTTGTAATTTGCCAGCATCAGCTTAAAACTCCTGTATGTTACATATCG TACCATTTAATCTCTTCACTTTCTTTGTACCTTCTGGCTGTATGCTTCTCTTTAACTTTTATAT TGTCATTTTATATTAATTTCTGTGTATATAATGTAAACCACAATTTTGAATAAAATTTAGTTCCTG CAGCTTGATTTAAATAGAGGACTTTTACTGGCACCTGCATCTCTCCAGATGCATGTACTCAGAGGGACT GTCAGACAAAAATAAATACATAAAAAACAAAGCAAGCAATGCACTATTAATTATACATTTTGATGTTCTA TCATTAATATTGTACGTACATTTTGGTTCACACAATTAATCCACTGTTTTCTCAGATGTAAAAATAA CCACATGCAGTTCTTGATTTA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_005786.6</u>
Summary:	This gene encodes a colon cancer antigen that was defined by serological analysis of recombinant cDNA expression libraries. The encoded protein is a member of the teashirt C2H2-type zinc-finger protein family and may be involved in transcriptional regulation of developmental processes. Mutations in this gene may be associated with congenital aural atresia syndrome. [provided by RefSeq, Jan 2012]
Locus ID:	10194
MW:	51