

Product datasheet for **SC213965**

ZNF280D (NM_017661) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF280D (NM_017661) Human 3' UTR Clone
Symbol:	ZNF280D
Synonyms:	SUHW4; ZNF634
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017661
Insert Size:	1352 bp



[View online »](#)

Insert Sequence:	>SC213965 3'UTR clone of NM_017661 The sequence shown below is from the reference sequence of NM_017661. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GACCTGGAAGACGAAAAAGAAAGAGTGAATTAGTCATTTTAAGTTTCAGTGTACCAACGATAAGGG CATTTGGAACAGTGCTATCAGGTGAGCTCAGTGGTGCTGTTGTAGGTCAGAAATGGAAATATGTAAGG GAGGTCACACATACACTTTACCTGTATGTTCAACCTATGTTATCAAAACAATCAATTCACCAATAATAG CATGATTAGTAGGGATCCCAAAAAGTTTTAAAAACACGAACAGGATTTTAATGATAATTAAATTTGC AGTGGAAGGTCTCATTTAATGGTTTTCAAGGAAATGGGATTTGGTTGCTGACATGAATTGATGATATT AGTAATATTTATAAAGCCTTTCAAACCTCCATCAATCCTAAGCTAAAAATCTTTATTACCTGTATATCC TTTTTCAGTTAACTGAGAGGAAGGGATTTGGAAACCATGTACTTTTGGGGAGTAATTGATTAACCAAT GGCTGATTGGCATTGTTAATGAAGGCTTTATTTGTGAGGATGATGCTGGTAAATGGAGCATGCTTAGAG TACTAAATTGATCTAATGAGAATTTGGATGAACATAAACTTAATTTGGATTTAATAACATTCAGT CAGACGCATGTAAACAGAATATTTGAATCTTTGTACCTCCATACAAGTGTTAGCTGCCAGGCTGTAAG CTTACCTTAATTAACCTTTTCAGTGAAAGTGAATTATTAAGATATAAATTTATATTTGTGCTTTTTGTC AGTGTGTAAGCTGTGTAGAAATCTTTGATGTATTAGTTGTATTAATGTAAAGTAGAAACCCATTGTTG AAACTCCTGTAGCTATTATGCTTTTAATATTGTTTTAATGATCTTCCTTAGAAATAGGCCCATAAAAAT GGTCTGGAAGCCAAACCAAGTATGGTATAATGTAGATATTGTAAAGCAGTAACTGAAACATGTCCT GGCATGTATTAGCCATGTTTAAGTGACTTTTCTGTAATTGTAAAATAAAAACTTCAAATGGGACCTAA AACAGTGATGTAAAGAACTGGTTTTGGAAATTTAGCCTAATTTATCTATAAGATGGCTGCTAAATTGA TTTTTCAGTTCTTTTTATCATCTAGAATATAATAGATATAGAAATGAATAATATGAAGAACAGTAGTTT GCTTTGAAATACTAATAAACTTTTATTTAAAATGCTTCATTTTTACTTCTTAAAATGTGCTTTGGATTCT TTAAATTTTGTTCCTGAATGTTCAATGTTTTAAATGGCGATTAATACTCTGCTGTATATAGTAGT TTTTGAGTAAATATTTGCAATAAAAAATCTGCCCCGAATAA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_017661.4</u>
Summary:	May function as a transcription factor.[UniProtKB/Swiss-Prot Function]
Locus ID:	54816
MW:	52.4