

Product datasheet for **SC213966**

ZNF280D (NM_001002843) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF280D (NM_001002843) Human 3' UTR Clone
Symbol:	ZNF280D
Synonyms:	SUHW4; ZNF634
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001002843
Insert Size:	1352 bp



[View online »](#)

Insert Sequence:	>SC213966 3'UTR clone of NM_001002843 The sequence shown below is from the reference sequence of NM_001002843. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GACCTGGAAGACGAAAAAGAAAGAGTGAATTAGTCATTTTAAGTTTCAGTGTACCAACGATAAGGG CATTTGGAACAGTGCTATCAGGTGAGCTCAGTGGTGCTGTTGTAGGTCAGAAATGGAAATATGTAAGG GAGGTCACACATACACTTTACCTGTATGTTCAACCTATGTTATCAAAACAAATCAATTCACCAATAATAG CATGATTAGTAGGGATCCCAAAAAGTTTTTAAAAACACGAACAGGATTTTAAATGATAATTAAATTTGC AGTGGAAGGTCTCATTTAATGGTTTTCAAGGAAATGGGATTTGGTTGCTGACATGAATTGATGATATT AGTAATATTTATAAAGCCTTTCAAACCTCCATCAATCCTAAGCTAAAAATCTTTATTACCTGTATATCC TTTTTCAGTTAACTGAGAGGAAGGGATTTGGAAACCATGTACTTTTGGGGAGTAATTGATTAACCAAT GGCTGATTGGCATTGTTAATGAAGGCTTTATTTGTGAGGATGATGCTGGTAAATGGAGCATGCTTAGAG TACTAAATTGATCTAATGAGAATTTGGATGAACATAAACTTAATTTGGATTTAATAACATTCCAGT CAGACGCATGTAAACAGAATATTTGAATCTTTGTACCTCCATACAAGTGTTAGCTGCCAGGCTGTAAAG CTTACCTTAATTAACCTTTTCAGTGAAAGTGAATTATTAAGATATAAATTTATATTTGTGCTTTTTGTC AGTGTGTAAGCTGTGTAGAAATCTTTGATGTATTAGTTGTATTAATGTAAAGTAGAAACCCATTGTTG AAACTCCTGTAGCTATTATGCTTTTAATATTGTTTTAATGATCTTCCTTAGAAATAGGCCCATAAAAAT GGTCTGGAAGCCAAACCAAGTATGGTATAATGTAGATATTGTAAAGCAGTAACTGAAACATGTCCT GGCATGTATTAGCCATGTTTAAGTGACTTTTCTGTAATTGTAAAATAAAAACTTCAAATGGGACCTAA AACAGTGATGTAAAGAACTGGTTTTGGAAATTTAGCCTAATTTATCTATAAGATGGCTGCTAAATTGA TTTTTCAGTTCTTTTTATCATCTAGAATATAATAGATATAGAAATGAATAATATGAAGAACAGTAGTTT GCTTTGAAATACTAATAAACTTTTATTTAAAATGCTTCATTTTTACTTCTTAAATGTGCTTTGGATTTC TTAAATTTTGTTCCTGAATGTTCAATGTTTTAAATGGCGATTAATACTCTGCTGTATATAGTAGT TTTTGAGTAAATATTTGCAATAAAAAATCTGCCCCGAATAA ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001002843.3
Summary:	May function as a transcription factor.[UniProtKB/Swiss-Prot Function]
Locus ID:	54816
MW:	52.4