

Product datasheet for **SC213967**

ELMO2 (NM_182764) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: ELMO2 (NM_182764) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: ELMO2
Synonyms: CED-12; Ced-12A; CED12; ELMO-2; VMPI
ACCN: NM_182764
Insert Size: 1340 bp
Insert Sequence: >SC213967 3'UTR clone of NM_182764

The sequence shown below is from the reference sequence of NM_182764. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGCTATGACTTTGTCTATCACTATGGCTGAGCCTGGAGCCAGAAACGACGGTACCCAGGAGAAGGGATT
TTGGGCCCAGGAGAAACACTTACATTCTGGTGCCTTGTCTTTTGGCTTGTACAGAATCTGTAGTGATTTT
GGTGGCCAGTAAATGCCAGCCATTTCTCAAACCCACCTCGGACCACCCAGAGTTTCTCTTGGTCCCTG
TCTACTAAGAGTCATGAAGGCAGGGTGTCTGCCACTCCATCACCATGAAGCCTGGGATTGGGCCACG
AGGAACAAACAGCAGATGCCCTTGCCTTCCAGTCCAAGAACTGCTTCTTGAATGGATTAAACAACAG
CCACTCACCTTTTCTCTGAGCCTGCTCTGATCAGCTGGATCCCCACGTGAGCAACAGCTGGCCCA
GGAAAGGCTGCCTGCAGAGGACAGGTGTGTTGGGCGTGTGAGAGCCTTGAAGTGACTACCTGTATCTT
AGATCTGAGTACAAGCCTGAGGCTTTTGTCTTTTGTCTTTTGTGATGAGGGCTCACTCCAGCTTCATATG
GTGCCAAGACGTTGCTGCTTCTGAGGTTGGCTTAACATCTCTGGTCTTTAGAGCCACCAGATCTCTCT
GGCCCATACAGATATCAGAGCAGACGGAATTTCTCCCTGCAAGCGCTCAGTCTCATCCAGCAAGTCA
AAGACCTCCTGGCCAAGTCTGCCTCTTAAGTCTCCAGGAACGCTGCAGGAAAACCCAGCTGAGGCC
TGGCCTAGACTGTGGTGAAGTCACTAGATTCTACTGCTTCCCCACATTAATACCTTTTCTTCTCT
CAGAGAGAAATCTCCCCTAACCTGAATTGCAGCCCTCCAGTTTGTCTTCTTGGCCTTCCAGACCC
CAGGAAGTTGGCCTTCCCTTCTAGTGCTATGGTTTCTGCCATTGGCCATGATTTAGGGAGCTGGCTG
AGGCCGGCTGAGGCCACCTGTGCCAGTGGGCTTCCCTGGTGTGCAGCACTTGTAAACCACACACA
CAGCCTCTCTCCCTGGACATACGTTAGCACATTGGCATTAGTATTGGTGGCCTGGCATGGTAGGTACT
ACCCAATGAAGAGTGTACTATATATTTTCACTACTATAGCCATACTTATACAGACGTGTATATATATT
TATATAAGATCTACCTATCTTAGGATGGAACCTGGGGAAAAATAAAATTGAGGGGAAGTAAAAAGTAT
GTAAACACTTCCAGTTGTGAGCCAAGATTGTAACCAAGAGAGCAGCCAGGAGCTTCTGTGAGTAACCATG
TTTTCAATAAATACTCTTTCATGTACAAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```



[View online >](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_182764.3
Summary:	The protein encoded by this gene interacts with the dedicator of cyto-kinesis 1 protein. Similarity to a C. elegans protein suggests that this protein may function in phagocytosis of apoptotic cells and in cell migration. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2015]
Locus ID:	63916
MW:	51