

Product datasheet for **SC213973**

C9orf150 (LURAP1L) (NM_203403) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C9orf150 (LURAP1L) (NM_203403) Human 3' UTR Clone
Symbol:	C9orf150
Synonyms:	bA3L8.2; C9orf150; HYST0841; LRAP35b
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_203403
Insert Size:	1330 bp



[View online »](#)

Insert Sequence: >SC213973 3'UTR clone of NM_203403
The sequence shown below is from the reference sequence of NM_203403. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TTGGATTCTGAATACTACTGCTTTGGCAGTGACAGTTTTTTGCATGGGACTGGTGTGAATGAAGTTG
TATTTATCCTTCTCTCCGCTGCTATATTTTGGTGTGATTTTTATTTTAATAAGATGACCTTTTTAAA
AGAAGCTGATTTTGAACTGCTTAATGGTATTGCTGTTGCTCCTAATACTTCTCATCTGAGCTGATTTA
TTTTCTCTGTTACATCTCTATTTTTTATTATTACAATGATTTTCTCCCTTCTTTTACAGTAGCACAA
ACAAAGTAGGGGAAAAGAATAAGCAATAATTATGTTTTGCTTTTGTTCAGAGCAATGGGTACAGGG
ATTACAAGAAAACTTTGCTAAATTTACAATAAACCAAGTCTGATAACAGTTAATGTTGCTGCTTGC
GTCCTTAAATGACTTAAGGTTTCATCTCCAGAAGATATTGAGATATATTACTGTTTGCATATCAGGTT
GAACAAAATGCTTGAACAAAATTAAGCATTATTTCTAGACCACAACATTCTAGTGATAATCAACCAT
ACTGTTAGTCACACACGCCACATGACAATCATCTGAAAGTCCCCAACTATACACATTCTGTTCAACCA
ACTCACCTTGTGTTGGTACTTTAGATTAACTTTTGTCTCAGTGATTCTAACTCAAATCATATGTTT
ATTTGTGAATAAAATGTATTATTCTCCCAATTTCAAGGAATAATTGAGATTAAGTATGGAGGCCATT
CAAGGCCAAAAATGCTCTAAATTTAATTGCATTTTGAGATTTTGTGGCATTACTTGTATGATTTTGT
GATTCTGATAGGATGTATAAATACTTAGATGCATACAAGCCGATGGATAGAAAGACAGATAAATATCCC
GATAGCCTACTATCAGCATCCTACTACTGCTTCCCTAAAGCAGAAATTGTTATTCATTTAAGAAGTCTAA
TATGATATACATTATCCTAGTATAAATAAATAATCTCTAAGTTTTCTGAAATGGAATATTTAACGCAC
ATAGACAGTGGATATTTACTTTAAATATAATGTTATTAATCCTCTGCATTTCTTTTATCTTGATAACTC
CTAGTGAGGGAATCTAATTAAGATTTTCATGCTGTTCTTTGTAATGATAAATGTTTCCGGACTAAACACT
CTATAATACAAATCAAAAAATCTTCCTTCTTATAACTGAGAAATGTTTACCCAAATTTCTTTGTTAT
TATATATTTACTACATGTTGAATTTATTATTGTTTATTTCCAAGGCATGTATTATTCAATTTCTAAT
TAAATATATCAGGGCTTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_203403.2

Locus ID: 286343

MW: 51.8