

Product datasheet for **SC213978**

LYSMD1 (NM_212551) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LYSMD1 (NM_212551) Human 3' UTR Clone
Symbol:	LYSMD1
Synonyms:	SB145
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_212551
Insert Size:	1164 bp
Insert Sequence:	>SC213978 3'UTR clone of NM_212551 The sequence shown below is from the reference sequence of NM_212551. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GACCAGGAGGATGAAATCTTCAAATCTTGATGTCCCCAGAACTGACTGAGAGAAGCAAGATGTTGAAAG
AGAAATTCAGGGGGAGAGGAGCCTGAGGTGAGGCTCAAGAACATGGCTTATGAGCCTGCCTCCCTCCA
CTCCAGCCATCTCCCCAGCCTCCTTGCTGCCAAAAATCCTTGGCCTGGGGCAGTTTTATTGGCAAG
ATTAGGGAATGGGAAGTGGACCTTTCCAGTTCTTGGGCTGAATCTAGAGTTGAGTCCTTTAGACTCAA
AAGGCTCTTTTTATATCCTTGTTGCCTTTCTATCCTTTTTACCTTTCCCTGGCCTTAGAACAGGGAAG
CAGAGACTCCCCCTCAGCTTCTTCTCATCTCCAACATGTGACATAATTTATTTGGTGCTAATATTGA
AGTAATATTTATTTGCTATATTTATTTCTTCCCACTCTTAGCTGAAAAATGGAATTTCTATAGCCCTG
GGGTGAGGAGAACCCGGGAACAGTGAGCCAGTGACTGCTGGTTCTCCACTTTCCCACTTTGTTAGCCA
AACATGTTATTGAGGATGAGGCAATCGGTGATATAGCCCAACACAGAGAATGGATGACTGTTGGCCTTT
TCTCAATCTCAAAGTGACAGATTTGGGCTGATTCTTTTTTTTTTTTTTTTTTTGGTACCCCATTCCTGTT
GAATGGGCTCAGTCTTGACAAACCAATTCCTGTAACCTTTTACAATCTTGACTCTGAGAAAAAGCA
GCACTGTAGGGACCTGAACCTCTCATTCTACCAGTTTACCCTTCTGACCTGTTCTGCTCAGCCTGTGGT
GAGGGGAACAGATGGAAGAAGGCTTTCCTTGGGGTTTGAAGAGGAAGCCAGTTGTGGGAGTCTCTTTGT
TACACCCCTCCCTCCAGACCTCTCTCTATGCCTTCTCCTCACTGCATATTTATTTGCGATTG
TAGCACTGATCTATGGTACAGTTGTTTTGTTTCTTAAACAGAAGTAGGTAAGGGAACCTCAGTAATAAT
GGAAATCTGCCTCAGTCTCCTGTTTGGAAATGATGAGAGGTGAGATAAGTAATGTTTATATACCTCCTT
TCTTCTGTGATGATTATAAGGTATACGAAGCTGGACAGATTAAACATTGAGCCCTGTGCA
ACGCGTAAGCGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```



[View online »](#)

Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_212551.5</u>
Locus ID:	388695
MW:	43.7