

Product datasheet for **SC213981**

CAMKV (NM_024046) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CAMKV (NM_024046) Human 3' UTR Clone
Symbol:	CAMKV
Synonyms:	1G5; VACAMKL
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024046
Insert Size:	1359 bp



[View online »](#)

Insert Sequence: >SC213981 3'UTR clone of NM_024046
The sequence shown below is from the reference sequence of NM_024046. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGGAGTCTCAAAGGGAGGAGGCCAGCTAGTAGGCAGCCTGGTGAGGGGGGCGAGGGGATGGGCAGGA
GGGTGGGAGAGTGGATGAGGGGCTTCTCACTGTACATAGAGTCACTGGCATGATGCCCTCGCTCCCCCA
TGCCCCACATCCCAGTGGGGCATAACTAGGGGTACGGGAGAGCAGTCTCGTCTCCTGTGTGTATGTG
TGTGAGTGGTGGCAGGCCAGTGGCAGGGCCGGCCCCAGCCCCTGCATGGATTCTTGTGGCTTTTCTG
TCTTTTGTAGCTTACCAGTTTCTGTTCTTGTGGGATGCTGCTCTAGGGATACTCAGGGGCTCTCTG
CTCTCTTCCCCTTCCCTTCTTGCCTCACCATTCCCCTAGGCAGGCCCTGCAGGTCCCACACTCTCCCA
GGCCCTAAACTTGGGCGCCTTGCCTGTAGAGCTGGTCTCCAGCGAGGCCCTGTCAGCGGTCTTAGGC
TCCTGCACATGAAGGTGTGTGCCTGTGGTGTGTGGCTGCTCTAGGAGCAGATACAGGCTGGTATAGAG
GATGCAGAAAGGTAGGGCAGTATGTTTAAAGTCCAGACTTGGCACATGGCTAGGGATACTGCTCACTAGC
TGTGGAGGTCTCAGGAGTGGAGAGAATGAGTAGGAGGGCAGAAGCTTCCATTTTGTCTTCTCTAAGA
CCCTGTTATTTGTGTTATTTCTGCCTTTCCGAGTCTGCACTGGGCTGCCCTGTACCCTGAACCTCAT
GAGCCTCTAAGGGAAAGGAGGAACAATTAGGACGTGGCAATGAGACCTGGCAGGGCAGAGTACAAGCCC
AGCACCCAGTGTCCCAGCCTTACTGGGTCTTACCCTGGGCCAAACAGGGAGGGCTGATACCTCCTTGC
TCTTCTAGATGCCCACCTCCTACAATCTCAGCCACAAGTCTCTCCACCCTAGGGGGCTTGTGCAT
GGCAATAACTCATAATCTGATTTGGAGGTTTGCCTTTACAGGGGCAGATTTTCTGCTCAGTTCAACAA
TGAAATGAAGAGGAACCTCCTCTTTCTACAGCTCACTTCTATCAGAGGCCAGGTGCCTCAGAGCCACA
TTGAGTTGCTTTTTCTGGGATGAGGAAGTAGGGTTAAACTCCCCAGTTTCTGAGGGAGGCTCCTGACA
GGTGCCCTTTGTGAGACCCTACCACAGCCTGGATAGGCAGCCACATTGGTCTCGCCCTTGTCTGGCAC
TCCGTGGTGGTCTGCCCTTCTCCCTGCATGCCTGTGGTCTGCTCTGGTGTGTGAAGGTCGGTGGGTT
AACTGTGTGCCTACTGAACCTGGCAAATAAACATACCCTGCAAAGCC
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_024046.5

Summary: Does not appear to have detectable kinase activity.[UniProtKB/Swiss-Prot Function]

Locus ID: 79012

MW: 49.3