

Product datasheet for **SC213987**

TNRC18 (NM_001080495) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TNRC18 (NM_001080495) Human 3' UTR Clone
Symbol:	TNRC18
Synonyms:	CAGL79; TNRC18A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001080495
Insert Size:	1325 bp



[View online »](#)

Insert Sequence:	>SC213987 3'UTR clone of NM_001080495 The sequence shown below is from the reference sequence of NM_001080495. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TCCACGGACGGCGTGCCCGTGCTCTGCTAGCCCCGCCGGCCCTGCGGGCCACCTGTGCCCCGAGGGCG GCCAGGGACCCATCTCCATCACTGCCATGGCGCGGAGACCACGTGCGTTGTGTGCATGCGAGCGCTCCT GCAGGCGTGTGCATGGGGCCAGGTGTGCACGTGGATGCCCGGGTGAGTGTGTGTGCATGCACGTGTGTG CACGCATGTGCACATGCCTCCAGCCCCACCTTCCAACCCCTCAGTGCCCCCAGGACAGGGGCCCCCTCTT AGCTATCAGGGTATGGCCGACCGGCCCTTCTGCCCAGCACGTTGCAAGCACTTGCCAGGCCGGCCCC TCCAGGCTGCTGCTGCGTGGGGGCCCGGGTGCCCCCAGGTCCATGCAGACTGGGGATTGGTGGGGAGG GGCGCTTCTAAGGAACCAAACTGACGCTCACTCTGGGCTTCCCAAGCACCCCTAGCATGGAGCCACCC CAAGGCCTCCACCGTGATGGGACAGGGCAGCCCGATGCCAGCTGGCCTACCTGACTGTGCCAGGGGC CCTGCCCCCACCTCTCAGGATGGCCTAGACTTGGGGAACAGAGCCGGGTGGGGTTGCAGCCCGGAGTG TCTGTCAAAGGCACCAAGGTGGAGAGGGGCCCGGCACAGGCCACCCCTGGTCCAAACCTCACACTACAGA AAACCCCAATGGTGCTGAAACTGTGCCCCGCCACGCCTGGCCCCCTCCCCACCCAGGAGGGAGGTGGCA CTTCTTAACCTGTACAGTTTTATTGTACCAAGAGACTCGCCCCGCCCTGTATCATAAGCCTTTAAATG GAGTCACTTTTTAATTATATATAAAGATAAATATATATAAATATATATAAACTTTTTAAACTG TGAAAAATAGCTATGAAATTATAAAAAAACATTCTGACGTGCAGAATATTATTTTTATTTCTGTGTT AGATTATAGTGTCTAGCACCGGCTTACCGGCCCTCCAGTCCCCAGCACACCCCCGCCACCCACCA AGTGTACTGTACTACCCCCCAGGATAGAGAAGTGTGTTAGGGAGAGAAGAGGGAGAGGCAGGAGCC GGCCCAAGCCCAGGGTCCTGCTTGGGCCCCAGAAAGCACTTAACCAGGCCCAAGCCTTCAAGGGAAA CCAAGGCCTCAACCAGACAATCTTGAGGGAAGGAAAAGCCAGACTTTGGGTTTGTGTTTGGGGGAATT ATTGGTTTTTTTTTTATGTTTCTTTTGAATTTGTTTGTGGCAAATCTGTGTGATCTTTTTTCAT AAAAAAAAAGAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001080495.3</u>
Locus ID:	84629
MW:	47.4