

Product datasheet for **SC214000**

ZNF182 (NM_006962) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF182 (NM_006962) Human 3' UTR Clone
Symbol:	ZNF182
Synonyms:	HHZ150; KOX14; Zfp182; ZNF21
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006962
Insert Size:	1346 bp



[View online »](#)

Insert Sequence:	>SC214000 3'UTR clone of NM_006962 The sequence shown below is from the reference sequence of NM_006962. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC ACTCGGAAGTCAAAGTTCATGGCACATTAGAGAGCTAATCAGCAGTATCTATTATGGACACTGAAAGAA AGTTGTGTCAATTTTACTCACATTTTAAAAAGTATATTTTGAAGTTTAAATATGGGAGTAAATTC AGCCTTTCTTCTTTTCTCAGTCATAACCCAAAAGCTACGGGAGAAAAAGCAGAAGTGTAAATCTCTG TATCTGACAAAATTAGGAGACAGCTGCAACCCTAATGAGAGGGCTGCCACAGGCAGAGGAAATGAAGCA CTGGTGCAGGAGATTTTCAAGTGTAGAGTGTCAAGGCTGCAGGTGCAAGACAGCACTGGTAAGGACTCT TAGCTGAGGTTCCGGATATACCCGAGGCGCGGAACCATAGATATAACAGTGGTAACAGTGCCTGGGCA GGTAGGAGGTAGGATGTTTCTAGACCATTGGGGTCTGAGGAGAAACAGGCTGGGGCTCGAGAAGGA GTCACCTCAGACAAGCCATATTTTGAAGTAGTCTACTGATGGGAGCAACAGTTATGGCATGCTGGGG GAAGAGAAGGGACTTTTCTTGTGAAGAGCCCTACAGCTGCCCTATTTACTGACTTGGGAGAAGTAAA GTCAGAAAGAACGAACAGCTGTGGGGGAAAAAATTCTCAGTTGTCAAAAAAAAAAAAACTGATCTGC TTAGAATGTGGCCCTGGTTTCTTTCCACATGAACCTCCAGCAAATAGCAAGTCCAGGAAAAAATCTT CTCATTCAAAGAGGAATTATTCAAAGTAATCAAAACCACAAATGCATTTGGTTTTTGGGCCAGCAATT CCACTTCTAGGAATGTACATACTCACACATGTACAAAATGATGAACATACAGTATTCATTAATGCAGCG TTTGTGAAAAACAAAGATTGGAAGCTGTTAAATGACACCATATCAGGGTCTGGTTCCCTAAATGTG GTATATTCATACAGCCTACTCTGGGTAGACAACAATGAATAAATAAGGAGAAGTTGGAACCATTTAT GTAGTGATGTGGGATAACCTTCAGTTTATGTTATTAGGTAGGAAAAATAAGGTTGAGGACAGTACTGC CATTTGTATATAAAAGGAGAAAAATGGTAAGTACATAATTGTGAATTCATAAGTTATCTCTGGAAAGA GACATGAGAACCTGGTAACCTGATTACCTTGGGAGGGGAAGGGGCAAGTGGTGGCTGGGGGATGGAAT TGAGAGGGAGACAAGTATCCTTTTATAGCTTTTGAAGTGTGGGTGTACTATGTACTATGTGGATGT ACTACCTTTTCAAAAATAAACATTTAAAAAGCTA ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_006962.2</u>
Summary:	Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This gene encodes a zinc finger protein, and belongs to the krueppel C2H2-type zinc-finger protein family. It may be involved in transcriptional regulation. Multiple alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, May 2010]
Locus ID:	7569
MW:	51.7