

Product datasheet for **SC214001**

HDHD2 (NM_032124) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	HDHD2 (NM_032124) Human 3' UTR Clone
Symbol:	HDHD2
Synonyms:	3110052N05Rik; HEL-S-301
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_032124
Insert Size:	1302 bp



[View online »](#)

Insert Sequence: >SC214001 3'UTR clone of NM_032124
The sequence shown below is from the reference sequence of NM_032124. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GTGGACCACATTCTGCAGCACCTATTGTGAAGCAGTGTGTGCATCTGAAGCAACTTGAAATGCAGCTTC
TTATTGTCTGGAATGAATCCCTTACCAACTCAGTGCCAGCATCGGTAGACACCAGTCAGTGCTGATCGC
TTTTTAACCTCTTTTGTGTGCATTAATTAGAAAAGAAAGGTATTGAATTGCGGTAGCCAGTAAGCCT
TGCTAATCTCTTTTATTTTGTAACTGAAGATGAGACCCAAAGAAAGGAAAGCTGAGATTTTGTGCCAT
TCCTTTTAAAATATTCATCAGGTTAGGTGGGGCTGTGGGGGAAAAGCTACTACAGGGAAGAGTGTCTC
TGCTGTCTCTTCACTGAAAAACAGGGAGGGGGGATTTCACTGTGAAGAAAGTTGAATGGTGGTTTTT
AAATTATAAAGTAATGTATTTAAAGGTGCATTAGGCTGTAGTTCTAATATTGAGTTCAACTGTGAAATC
CATCAGATGTGCCAAATGGAGAAGACAGAAAGCAACAAAGTGAATTGTTCTTTAGCCCAAGTGGTACAG
TGAATTTGCTTTAACAGATGTTGAAAACTAAATTTTCTACTGTATTCCCAGCACGGGTGACTTCTTTTT
CTTTCATTAGCCAGAGATGACTAATTTAAATTTAGAACCAGATTTTAATTTAAATTAATTTCCATT
AATAACCTATTTCATTGCAGATACCTATTATACTGTGTAACAGTTGTTTTGGAAATTTTATGTAAAATTA
AACTATCAGTATTTTACAGATGTTTTAATTAGACATTGTTATTAACAGGAACAGTGCAGAACTAGAA
TCAAGCCTTATAATATCTTATAGACCATGCATTTTTGAAGTTAGTGTCCACTAGGGTCCCTATTAAGTGT
ACATTTGCAAGATTTTATTATTTTTGCCTCTGACACTATGGGAAAAATTTTTAGAAAGCTATTGGGACA
GATTCAGCTTTTATGCACTTGGTTACTACAGCTGTAAGTAAATGAAATCTCGTCTTTAGCATGGATTATT
CTTCTCATGTTAAACCCACCAAAAATAAAGGGGACTAAATAGGTAATGATTTTCTAGTGCATTTGCATA
CTGTGATAATCCTGGCCCTTGCAATAGTTCTACAGGGCTCTTGGGCATTGAATTATTAGGATGTAATTG
TACATCATTGTAGTGTTCACCTTATTGAAGCTCACTCTGATGTTAATGAGCTTCGGGTTTTGATGCTTG
TTTAGAGATCAGCAGTCTTGGATGGGAGGGAACAAAGCTAAATAAATGTTAGTTTGGTGA
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_032124.5](#)

Locus ID: 84064

MW: 49.6