

Product datasheet for **SC214018**

MCF2L2 (NM_015078) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MCF2L2 (NM_015078) Human 3' UTR Clone
Symbol:	MCF2L2
Synonyms:	ARHGEF22
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015078
Insert Size:	1369 bp



[View online »](#)

Insert Sequence:	>SC214018 3'UTR clone of NM_015078 The sequence shown below is from the reference sequence of NM_015078. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CGCCCGAGGACCTCCGCCAGGAGAGCTGACCTCCCTGCGGACGCCCGCTCCTCGGCTCCAGAGCGC CCGCATTCCCGGAGAGGCGGTGTGGGGGCCCGGCCCTGCCAGCTACGCAGAAAGCAGCCGGAGCCT CGGCGGCGGCAGAAAGGGACAACCAGGGCCTCTCCGAGGAGCCGAGGGGTGTCTGGGTGCGCGCC TAGCTCCGCACGGGGACCTCGGAGCTGCTCTAAGGCGCTGCAGAGGCGAGCAGAGCCCGAGCCAC GCCTTCTCGACCGCGCACTTCGATATTCGGAGCCGGGCAATTCTTTGCTGCGTGGGCTTCTCTGTGCTC CCACGGTAGGATGTTTAGTAGCACCCCTGGCCTTACCCACTAGGTGCCAGGAATGCGCCACCCCATCC TCCACCCCGCCCCAAATCGTGACAATCAAAAATGCCTGCAGACACGCCCGCATTTCTCCAGGGCGGGG TGGGAATCGGTTGAGAGCCGCTTAGCCCGAGCCTTGAGAGGAGCCGAGCCGCTCAACCCGCGGGGGC CGCAGACTGGGAGCTCCCGGTCCGCCTCCAGCATCCCTGCGAGCGTTTATGGGGTGTTCGTGTTAGTG CCAAGATTGCTTCGTTGTAGAGAGAGTTCGTTCCAAGTTACTTTCTGAGGTATTTTATGTATCGATTTA TTAGTTTTTAAATGAGTTTTGTTAGTTTCAGTTGTATTTTATTTTATTTTATTTAGTTTATTAGTTGATTTTATT TCTTTGTAGCTTTCTGGTTATTTAATTTTATAGTTTCACTTACTGGTTTATAGTTACTTTATTTTCAAA GTTATTTAGTTGTTTATTTTCACTTATTTATATAGTTGTTTGTGTTTATAGGAGTTACAGATGTTCAAA TTAATTTGCTTGGAAATTTATTTATTTATTTATTTATTTTTCGAGACGGAGTCTCGCTCTGCGCC CAGGCTGGAGTGCAGTGGTGTGATCTCGGCTCACTGCAAACCGCTCCCGGGTTACGCCATTCTCCTG CCTCAGCCTCCTGAGTAGCTGGGACTACAGGCGCCCGCCACCACGCCCGGCTAATTTTTTATTTGGATT TTTAGTAGAGACGGGGTTTCCCGTGTTAGCCAGGATGGTCTCGATCTCCTGACCTCGTGATCCCCCTG CCTCGGCTCCCAAAGTCTGGGATTACAGGCGTGAGCCACCGCGCCAGCCGGAATTTATTTTAAATT ATCTTTACAATTATTATCTGAGTTATTTACCTTCATAGTTATTTCACTTTATCTTATTGGAGTTTTTGA GTTATTTATTTTACAGTAACCTATTTGACTATTAACACTCACTGGAAGTTCATGAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_015078.4</u>
Summary:	Probably functions as a guanine nucleotide exchange factor.[UniProtKB/Swiss-Prot Function]
Locus ID:	23101
MW:	51.5