

## Product datasheet for **SC214034**

### **LPPR2 (PLPPR2) (NM\_022737) Human 3' UTR Clone**

#### **Product data:**

|                           |                                               |
|---------------------------|-----------------------------------------------|
| Product Type:             | 3' UTR Clones                                 |
| Product Name:             | LPPR2 (PLPPR2) (NM_022737) Human 3' UTR Clone |
| Symbol:                   | LPPR2                                         |
| Synonyms:                 | LPPR2; PRG4                                   |
| Mammalian Cell Selection: | Neomycin                                      |
| Vector:                   | pMirTarget (PS100062)                         |
| ACCN:                     | NM_022737                                     |
| Insert Size:              | 1349 bp                                       |



[View online »](#)

**Insert Sequence:** >SC214034 3'UTR clone of NM\_022737  
The sequence shown below is from the reference sequence of NM\_022737. The complete sequence of this clone may contain minor differences, such as SNPs.  
**Blue**=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CGCAGAACTGCGCCCGCGTGGCCACCTGATCCCCAGCTGTGTCTCTCCAGGGCCCCAGCCATGTGTT
CGTCGCCCCGTGTGCCCGTCTCGATTGAGGTCTGAGCCGACGCCCTTGCCCTGCCCTACCCCTGC
CAGCGCCACCCCGAGCCAGGGCCCTCGCCTTCTCCCTGGACCTGGGGGGCCAGGCGGGGGTGGTG
GACGTGGCCGGAAGCTGCTGCTGCCACGCCCTGCTGCGGGACCTGTACACCCTGAGTGGACTCTATC
CCTCCCCCTTCCACGGGACAATTACAGCCCTTACCTGTTTGCCAGCCGTGACCACCTGCTGTGAGGCC
CGACCACCCAGCAATCTGCCAGTCCCACTTCTTCCCTGCCACGCGTGTGTGTGCGTGTGCCACG
TGAGTGCCAAAGTCCCTGCCCCCAAGCCAGCCAGACCCAGACATTAGAAGATGGCTAGAAGGACATT
TAGGAGACATCTGCCTCTCTGGCCCTCTGAGATATCCCGATGGGCACAAATGGAAGGTGCGCACTTGCC
CCTACTATTGCCCTTTTAAAGGGCAAAGCTTGACCCCATTTGGCCATTGCCTGGCTAATGAGAACCCCTG
GTTCTCAGAATTTTAACCAAAAGGAGTTGGCTCCAACCAATGGGAGCCTTCCCCTCACTTCTTAGAATC
CTCCTGCAAGAGGGCAACTCCAGCCAGTGTTCCAGCGACTGAACAGCCAATAGGAGCCCTTGGTTTCCAG
AATTTCTAGAGTGGGTGGGCATGATTCCAGTCAATGGGGGACCGCCCGTGTCTAAGCATGTGCAAAGGA
GAGGAGGGAGATGAGGTCAATTGTTGTCTTGTGAGTCTTCTCTCAGAATCAGCGAGCCAGCTGTAGGGT
GGGGGGCAGGCTCCCCATGGCAGGGTCTTGGGGTACCCCTTTTCTCTCAGCCCTCCCTGTGTGCG
GCCTCTCCACCTCTCACCACTCTCTCCTAATCCCTACTTAAGTAGGGCTTGCCCACTTCAGAGGTT
TTGGGGTTCAGGGTGCTGTGTCTCCCTTGCCCTGTGCCAGGTCTATCCAAACCTTCTGTATTATTATT
AGGGCTGTGGGAAGGGTTTTTCTTTTCTTTTCTTGAACCTGCCCTGTTCTTCACTGCCCCCATGTC
CTCAGCCTCATACAGATGTGCCATCATGGGGGGCATGGGTGGAGCAGAGGGGCTCCCTCACCCCGGCA
GGCAAAGGCAGTGGGTAGAGGAGGCACTGCCCCCTTTCTGCCCTCTCATCTTTAATAAAGACCT
GGCTTCTCATCTTTAATAAAGACCTGTTTGTAAACAGAA
AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCC
CAACCTGCCATCACGAGATTTGATTCCACCGCCGC

```

**Restriction Sites:** SgfI-RsrII

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_022737.3](#)

**Locus ID:** 64748

**MW:** 47.7