

Product datasheet for **SC214064**

JRKL (NM_003772) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	JRKL (NM_003772) Human 3' UTR Clone
Symbol:	JRKL
Synonyms:	HBMJG
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003772
Insert Size:	1367 bp



[View online »](#)

Insert Sequence: >SC214064 3'UTR clone of NM_003772
The sequence shown below is from the reference sequence of NM_003772. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAACAGAAGATGACAAAGTCAAGTCAATTAATGTCAATTTTATTGTTCTGCTCATTGTGTTTGTG
ACAAACTCTTTGCAATATGGCTTAATTTTCTTTGTGTTCTGAATTCTCAGACTTGGTCCTGTGAAATAC
AGGCACAAAATGTATCTGAAGTGGTTTGAGGATTATGTGTTTTCATCATCTGTGCTTTTGTCTTTTA
TTTGTACAGATAATCAGAAGATGATACTGAATAGATATAAATTACATGTACACATGTATTCACTTTTA
GAATCTGCAATTATACCTTCTGTAACTGGCATTCCCTTAATTTTCTAGTGAAAGTTAGAGATACTG
AACAGACTGAAGCACTTTTCTGAAATCTTTGCTTGATTTATGAAGGCTGCCATAGTTATCTTTTCTTG
TGTTAACCATCTTAAATGATGTTTGTATATTTAATAGACTGATAGGATGAGAAAGATTTATATTATT
AGATTTCAAGGATGATTTATAATAATTCAAAAATGAAATCAATAATGGGGAATAATTATGAATTATAG
AAATTATGCCTTCATTCTCTTACATTTGTGTGGTTGCAAGAGGGGAGATTATTCTGGAATGAAGTA
AATATGGGAAGTTATTGCCAAATGAGAGAGAACTATGGGAAAGCTGATCTATAAAGAGGCATTCTGAT
CAATTCATTTGTAGGAACTGGGAAATAAAACCTGGGGAACCTTAGGTTATTTATACAAAGGGAATAA
ATAGGCTGATTTTAATTTGGTAAGTTGATCTTTTATTATGAATTTGGTAATAGTATAGGTTTATTATT
TATTCATCTAATTTATAGTACAGGTTTGTAAATGTTACATGTGATGATGAGCTCCACCTTATATG
GGGGAACATCTTGGGAATTTGAGATTTAATAAGTTTTTTTTTTTTTTTTTTAGTGTTTTACTGCATA
CTCACAAATGTTGTCTATAATTTGAAAAATATTGTCATATCTGGCCCTTTGATGAGAAAAGGAAATTAC
AATAATAAAGTTTATGATTTTAAATAAGTCATATGTTTGTATCTGTTTATGAAGAAAGCAGAAATA
ATTACTGAAAGTGCCAGACACTAAGGAATATTATTGTTTATTATTTAATACATATAAAAGGGATTAA
TCTGCTAAAATGTAATCTAAATCAGAATTTTGATAAATTTTTTTGTAACTAAGTATGTTTATTCAAG
ACATTGAACTACTTTGCACATATGAATATTAATGTAAGTTGTAATTTAAAGTAAAGTGTTCATGC
TATTTCTGTTTGGCCAAAAATTTTAAAAAATAAATTACAATTGTTCTCTATTA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_003772.4](#)

Summary: The function of this gene has not yet been defined, however, the encoded protein shares similarity with the human (41% identical) and mouse (34% identical) jerky gene products. This protein may act as a nuclear regulatory protein. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, May 2012]

Locus ID: 8690

MW: 54.6