

Product datasheet for **SC214078**

FUBP3 (NM_003934) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FUBP3 (NM_003934) Human 3' UTR Clone
Symbol:	FUBP3
Synonyms:	FBP3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003934
Insert Size:	1360 bp



[View online »](#)

Insert Sequence:	>SC214078 3'UTR clone of NM_003934 The sequence shown below is from the reference sequence of NM_003934. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CAGGCGCAGGCCACAGCCAGGAGCAGTAGGACAGCGTCCTCGTGGCCAACCTCTCCCTCAAAATCATT GTCAAAGAAAACCTATTTGTAACAGAGGTTTTTACATTTGCAAACCTTTTGATGAAGAACTGTTGTTTT GTTCTGTGAAACACTATATTTAGATTAACGGAATTTTTCTAAAAATCAGGAAAATTAGTTACTAAAAAT TGCTGATCATTTTTGTTTCATTATTTTTGTTATTTCAAATGTATAAGCTCTGGGATTCTTTTTGGAGCA ATACCTACAAAGTCAGGCACCAGAATGTGCCTCAGAGCTGTGACATTTCAACATGATGGTTTTGTTTG GTTTGTTTTGTGTCCTTTTTATTTGCAGTTTTTCTGTTGCAACAGAAAGTGGCTTGAAGTCTTAGG TGGTATGTAACAAATCTTTTTAAAAATTTAAAGCAGTATTTAAGTATTCCTAAATGTGTAATTCAT TTAATGTTTTACTTCTAATTTCTTGATCTTGGCTGTCTGGTTTTATTGCATTTTTAAAAAACTGAAC TATTATGTATTGTAATGAATGATGTGAATTCTGAACCTGAGACAGATAATTGTATTGCTGTCCAACAGGG GTTGGGAGGGGCATTTTGTAGGGTTTGTATAATTTCTAGAGTTTCAAAGGGGTAATATAAAAAATGTGCT CGTGTGTTTAGCTACGTTTTTCTCATGCCTCCATTACCAGTTGCAATTGTGTATCTAAGACCTCCAAG CTTGTTTTTAAAAACAAAACAAAACCTTTCTCTATTATCTCAGAAATATAAATACTGTTATTTTTAAT ATTAAGAAATTCATATAACTTTTAACACACACTGTGGAACATTAAACCGTATAAAAAATGTAGTTATAA CTATTTTTAATCCAAGTGTTGTTTGCCTTTTTCTTTTAAATATTTTCTTAATATTTTGTGAGATTAA CTAGATGAATAAATAAATCTAGTATTAACCGCATTATGAATTAATAATTTTGATTAATGAAAGGGAT AATATGATTTCCAGTGTTTACTGTAGTGTATCTTGACAGATAACATGTATTTTTAAAGGAAAAAAG GGAATTGAAGCTATTTTTTCTTGCAATTTCTAATTGACCTGAGGGACATTCCGTTTGAAATGTACTGAAG TTACAGTTTCTGGTTTTTCTCCTTATTTTCTTATAATGCTTGAATGTCTAACTATTAAAAAAGACA ATTGAAAAATGTTATGCATGGGGTTTTTAAGAAAACAAAGTGTCTTTTTATTGACTGACAATTCATT TTACACTCTATATAATAAAATCTCCACAAGGCATCTTGTGGGCAAAGTC ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_003934.2</u>
Summary:	May interact with single-stranded DNA from the far-upstream element (FUSE). May activate gene expression.[UniProtKB/Swiss-Prot Function]
Locus ID:	8939
MW:	52.5