

Product datasheet for **SC214081**

LG13 (NM_139278) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	LG13 (NM_139278) Human 3' UTR Clone
Symbol:	LG13
Synonyms:	LGIL4
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_139278
Insert Size:	1360 bp



[View online »](#)

Insert Sequence:	>SC214081 3'UTR clone of NM_139278 The sequence shown below is from the reference sequence of NM_139278. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAACGCGATCGCC AGACACATTGTGGTGGATCTCAGTGCCAGGGGGTGCCGAGGCCTCTGGTCTCCTCAGGGTGCCCACTG GAGGATGGGTGGGGGCTCTGTGAGCACACGTATGCCTGTGTGAAGACAAGCGACCAACCTGCATTCA TGCCCTTGATCCGCACACGCCAGTGTGTACAGGCGCACCCCGTGGACTGGCCAGGAGGACGTGCAT CATCACTGCAATCAGCATGAACACCAGGACCCCTGGGACACCTAACGCCAGTTCCTGACCCCAAC TGCAGGCCCCCTACTCTGCAGCCTTTCCCATGCTCTGCCCCCTGTCTCATGCACGAGTGGTGGGGG TGCTGGGAGCGAGGGGGCGCATCCTCTGTCTTGTCTTCTCTTCTGCCCCCTGGTCTTTCCTG TTGGTGTCTCAGGCATTGTTACCTGCTCTTGTCTCCCGCTGTAGCCACAGCCATATCTGTCATGCTT CCCCGGGGCCACTCACCCTGCCCCACCCCGTCTGCATGCTCACATGGACGCGGACGACACACTCATT ACACACTCACACGCCCGACTGCACCCGCTGTCCATCTATAAACATGTACAGGACGCTGTGCATACAAG CAGATTACGACTGCCTACCAGGCTCTCTGTTACCTTGCTCTCTGCTTGGAGGGGGCCCCCTCTGCT CACCCCCAGCGTGCCCCCCCCGGAAGTATAAGGATGAGAATGGTGGTGCCAGCTTCTGAGGGCCTGCT TGGCCTCTGGGGGCAAAGCCCCCTGCCCCGAAGCAATAACAACAGCAGCAGAAGCAATCCCGGGGAGC TGCCAGGGGTGTCTGTCTCTGTTCTCTTGGCGTTTCCTTCATGGCCGTGAATCCAGCTCACCCATGG ACAGAAGACCTCCCTCATCTGAAAGCCTTTCTGCCTTCTTCTGTAGCTGTGCTTTCTCAGGCTTG GACTGTCCCCTGGTTCACCCTGTCTTTGGGATGTGGCCCTTGAGCTGGGGGGCTCCTTGGTGGATTA GAGAGCTGGCCGAGCTGCAAGCACCTGAGTTACCTGCGAGGGTGCCAGAGTGGAGCTTGGCCTGAGGG CGTGTTGCAGCCCTGGCCCCACCAAGTCCACAGAGGCCCCAGGAGCACTGCCTTTCAAGTGACCGGGTG GGCCTGGTGGGAGGGGGCCCACTAGCTGAGCTTCTGCTGCCCTCTGCTGTCTGAGATGCCAGGCT AGTGGAGAGGGCACACCATGTGCCCAAATGGTGTGATGCTCCTTGTAGCAGGATACCAGGGTACTGGGG GGCAATGCTATGATTAACCTGCTTCAAATAAAAAGTTCCCGCCCCATAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_139278.4</u>
Summary:	May participate in the regulation of neuronal exocytosis.[UniProtKB/Swiss-Prot Function]
Locus ID:	203190
MW:	47.9