

Product datasheet for **SC214083**

CD16 (FCGR3A) (NM_000569) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: CD16 (FCGR3A) (NM_000569) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: FCGR3A
Synonyms: CD16; CD16A; FCG3; FCGR3; FCGR11; FCR-10; FCR11; FCR11A; IGFR3; IMD20
ACCN: NM_000569
Insert Size: 1283 bp
Insert Sequence: >SC214083 3'UTR clone of NM_000569
The sequence shown below is from the reference sequence of NM_000569. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAATGGAGAAAGGACCCCAAGACAAATGCCCCATCCCATGGGGTAATAAGAGCAGTAGCAGCAGC
ATCTCTGAACATTTCTCTGGATTTGCAACCCCATCATCCTCAGGCCTCTCTACAAGCAGCAGGAAACAT
AGAACTCAGAGCCAGATCCCTTATCCAACCTCTCGACTTTTCTTGGTCTCCAGTGAAGGAAAAGCCC
ATGATCTTCAAGCAGGGAAGCCCAAGTGTAGCTGCATTCTAGAAATTGAAGTTTCAGAGCTACACA
AACACTTTTTCTGTCCCAACCGTCCCTCACAGCAAAGCAACAATACAGGCTAGGGATGGTAATCCTTT
AAACATACAAAAATTGCTCGTGTATAAATTACCCAGTTTAGAGGGGAAAAAAAAACAATTATTCCTAA
ATAAATGGATAAGTAGAATTAATGGTTGAGGCAGGACCATACAGAGTGTGGAACTGCTGGGGATCTAG
GGAATTCAGTGGGACCAATGAAAGCATGGCTGAGAAATAGCAGGTAGTCCAGGATAGTCTAAGGGAGGT
GTTCCCATCTGAGCCAGAGATAAGGGTGTCTTCTAGAACATTAGCCGTAGTGAATTAACAGGAAAT
CATGAGGGTGACGTAGAATTGAGTCTTCCAGGGGACTCTATCAGAAGTGGACCATCTCCAAGTATATAA
CGATGAGTCCCTTAATGCTAGGAGTAGAAAATGGTCTTAGGAAGGGGACTGAGGATTGCGGTGGGGGG
TGGGGTGGAAAAGAAAGTACAGAACAAACCTGTGTCACTGTCCCAAGTTGCTAAGTGAACAGAAGTAT
CTCAGCATCAGAAATGAGAAAGCCTGAGAAGAAAGAACCAACCACACAGGAAAGGAAAGCGCAGG
AGGTGAAAATGCTTTCTTGGCCAGGGTAGTAAGAATTAGAGGTTAATGCAGGGACTGTAAAACACCTT
TTCTGCTTCAATATCTAATTCCTGTGTAGCTTTGTTTCATTGCATTTATTAACAATGTTGTATAACCA
ATACTAAATGACTACTGAGCTTCGCTGAGTTAAGTTATGAACTTTCAAATCCTCATCATGTCAGTT
CCAATGAGGTGGGGATGGAGAAGACAATTGTTGCTTATGAAAGAAAGCTTTAGCTGTCTCTGTTTTGTA
AGCTTTAAGCGCAACATTTCTTGGTTCCAATAAAGCATTTTACAAGATCTTGCATGCTACTCTTAGATA
GAAGATGGGAAAACCATGGTAATAAAAATGAATGATAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI



[View online >](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_000569.8
Summary:	This gene encodes a receptor for the Fc portion of immunoglobulin G, and it is involved in the removal of antigen-antibody complexes from the circulation, as well as other responses, including antibody dependent cellular mediated cytotoxicity and antibody dependent enhancement of virus infections. This gene (FCGR3A) is highly similar to another nearby gene (FCGR3B) located on chromosome 1. The receptor encoded by this gene is expressed on natural killer (NK) cells as an integral membrane glycoprotein anchored through a transmembrane peptide, whereas FCGR3B is expressed on polymorphonuclear neutrophils (PMN) where the receptor is anchored through a phosphatidylinositol (PI) linkage. Mutations in this gene are associated with immunodeficiency 20, and have been linked to susceptibility to recurrent viral infections, susceptibility to systemic lupus erythematosus, and alloimmune neonatal neutropenia. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2020]
Locus ID:	2214
MW:	47.8