

## Product datasheet for **SC214086**

### CHCHD7 (NM\_001011667) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	CHCHD7 (NM_001011667) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CHCHD7
Synonyms:	COX23
ACCN:	NM_001011667
Insert Size:	1535 bp



[View online »](#)

**Insert Sequence:** >SC214086 3'UTR clone of NM\_001011667  
 The sequence shown below is from the reference sequence of NM\_001011667. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACAGGGAAAGGTGTTCCACTTACTTCTTGAGGTACAAAACTGCCGAGATTCTGGAATTCTATCGTGA
TGCAGAGAAGAAAGAACGGAGTGAAGCCATTTATGCCTACGGCAGCAGAAAGAGATGAAATCTTGAGAG
CAGTGGGAAATATGCCCTATTGAATGTTTGCATTAATAAGTGTATATAACTTAGAAGCAGATGAATAT
TTCTAATAAATGATTGCTGTAATATTTAAGACTGTACACCCCTCACCCAGACAGACCTTAAGTTCTTCA
AGTGGAGACAGTGAAGTCAACCCGTGTCCTTTTGTCTGCTCAGTGCCATGCCGATGGTATGTTGCT
GTTGGCTGTGTTGGCATGAGTTGCATGACTTCTGGAGGCATGGAGTTAGGTAAGGCTACATGAGA
AATTGAGCTTTTCCACTGGGTTTTGAAAGAAGAGTATGATGTGAACAAGTAAAGACTGAATGGGGCTGA
GATGAAGGCAATGTTTCCAAGGAAAGGAAATGTTATGAGCAAGAGTGTGAGGCAAGAGAAGCTGGAACC
ACATTCAGAGAGTATCCTGTAGATTGCTCCACCTAGAATCTCAGGTGGGTGGAGCAGTGGTGGGAGAAG
ACTGGAAGGTAAGTTGAAGTAAGGAATGTGTGGTGGGCCTCAGATCCCAGGCTCATTCCCTCAAATCA
CTTCTTACTTCCCTCACTTATCTTTGTTTAAATAAGGTTAGTACATTCAGTGGGCAAATGGGTTTTT
CTAAATAAATGACATAAAAAAGAACTGGTTGTCTTTCTTAAAGAAAGTCTTAGTTACAGGATTCTGTG
ATTAATAATCCTTAAATATTGGGTTGCCTTCTGCAGCTGTAGTATCAGTTTTTGAACCACAGTAATGGGA
AGAAAGACAAATGGATTCCCTTAGAAGTAAGATTTCTATTTGCAGGATGAGTTGGGCAGGGAAAAGGGT
CAGGGTTCATCAGGTGAAGTCAACTGGGATGAGACTAGAAGTCACTTTATGATATAAACACAATAC
GATTATTCAATGTGGTACTGGGTGACTGTGAAGCAGCCCTCCTTATGCCCACTGCCTTTTGAATC
GTTTGTGTTTTATTGATGGTAGTTTTATGAAGCATACTATTATTGAATGAAATTTAATGTGACTTAAAA
ACATTGCTTTTGTCCCTTCTTTCATCTGGTCTTGGGTCAAGAACATTGTTTTAATGGCTGCCGACAAT
GAACTGTCTGTCTGAGTCTAAAACCAAGCTCAGGTTTCTAAGCCACATGACCTTGATTGTTAAATGTAG
TTATATTTGAATAAAAAATGAATCATTTCAAAAATTTCTTTGATTTGATTGACTGCTATAAAAAATGAT
GGATTAATAATAGGACCAGAAAAGAAATGAATGCCTACCTGTGTAAGAAGTATTTTTGTATTTTAAGCAC
ACTTGTAAAGGGCTTGGAAATTTAATGTCACCTCTTGAAGGTACATACCAATATCAACATTCATCTT
GTAAGTGTGCTCTGA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_001011667.3](#)

**Locus ID:** 79145

**MW:** 59