

Product datasheet for **SC214090**

MORC3 (NM_015358) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MORC3 (NM_015358) Human 3' UTR Clone
Symbol:	MORC3
Synonyms:	NXP2; ZCW5; ZCWCC3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015358
Insert Size:	1366 bp



[View online »](#)

Insert Sequence: >SC214090 3'UTR clone of NM_015358
The sequence shown below is from the reference sequence of NM_015358. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GAACAAATGAGTGAAATCAGTAGTACTAAAGTATATGTTATGTAAGATAAAATATTTGCTCAATTCTT
TTGTTGTACAGCTTTCAAATATAAATTAATTTTGTATAGATATGATAGGCAACAGACTGAAAACC
ATAATCTTTACTGTATTCTATGCATTCAAATGTGGTCACAAATATTGTGGACACATTATCTTATGTTTT
GAAATACCTGTGAATTGTTGGCATTGAGCAGCTGAAGCTAACTCATGACTCTGTTTTGAATGTAAATAT
TTGTAATTAAGCCTGCACATATTTTTTATTGCCCTAGAGTACTCAAGTGTTTTACCACAGAGCTTTT
CAGGTTGCCCTAAGCTTTGTGCAATTTTTCTGGTCCCAAAGTGATTTTTCTCTAAGTCGAGGGC
TATGCCATAATACAAATGGAAATGTTACCTTTGATTTTCTATAAAGGAGTTAAATAGGATTTTTAAA
TAATGTAGTAACACTCCTGATACAACTCTGGTTATAAGTGAATTGAGCATTAAATGTTCTTTGTATAA
ATTCTGTCTGAAATATTTTATTCATGAAAATAAGGTAAAGCAAAAACCAACTCCATTTTGCCAGGATT
TTGTTGTGCTGAGATTGCCAATCACGGTTAAACACAAAGTGTTTTATAGAATGAAGAAACAGTCTTTA
ACAGAAAAAAGGTATTGAAATATTAATGGCCACCAAGTTTACTTTGAAGCCCATTTTTGGCTGTTTAG
TCAGCATGAAGTGGGCATGATAATTTTTTAATATTTCTTTTGTGAAATTTCTGTACAGCCTTTTGTA
GGATTACTACAGTTAATGTGAGTTGAGGAAGACAGTCTTTCTCAAACAATACTACATACCTTTTATTA
TATTACAAACCTAAGTGTTTTATCCTAGAGTCAAGGAACAAATTTCTTAGGATTATAGTATTACAT
GCCATAAAATACTATGCTTTATTGGTCCCATGTTTGTGCAATTTTAAAGAGATGGCTTTCTATTAAGT
ATAAACTATGTATATATAAGAACCATATTTCCACAATAAAATGTGCATTATTTTTCCAAAGTTTTA
TCATTGCTATTTATTTTACCTTTGTTTTTGAACATTCAATCCGTTTCATTTTGTATGTATGCTTAATAC
GTGTCGGTCATATACAGTATTGAATTTTTACTGTATAGTAATTCTGGAAGAGCAAAATAAATGAAGATT
GTTTTATTTGCCTGATAAAGTAATTGAAAGTGTATTTTTGTTATGAAGCTGTTTTCTGTACAAATTG
TAATTTCCCAAATTTTAAATATCTTATAATAAAATAAAATATATGATGGCTAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_015358.3](#)

Summary: This gene encodes a protein that localizes to the nuclear matrix and forms nuclear bodies via an ATP-dependent mechanism. The protein is predicted to have coiled-coil and zinc finger domains and has RNA binding activity. Alternative splicing produces multiple transcript variants encoding distinct isoforms. [provided by RefSeq, Feb 2016]

Locus ID: 23515

MW: 53.1