

## Product datasheet for **SC214107**

### **FHL1 (NM\_001159704) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	FHL1 (NM_001159704) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	FHL1
Synonyms:	FCMSU; FHL-1; FHL1A; FHL1B; FLH1A; KYOT; RBMX1A; RBMX1B; SLIM; SLIM-1; SLIM1; SLIMMER; XMPMA
ACCN:	NM_001159704
Insert Size:	1382 bp



[View online »](#)

<b>Insert Sequence:</b>	<p>&gt;SC214107 3'UTR clone of NM_001159704 The sequence shown below is from the reference sequence of NM_001159704. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TATTGTCCCGACTGTGCCAAAAAGCTGTAAACTGACAGGGCTCCTGTCCTGTAATAATGGCATTGTAAT CTCGTTCTTTGTGTCCTTACTTTCTGCCCTATACCATCAATAGGGGAAGAGTGGTCCTTCCCTTCTTTA AAGTTCTCCTTCCGTCTTTTCTCCATTTTACAGTATTACTCAAATAAGGGCACACAGTGATCATATTA GCATTTAGCAAAAAGCAACCCTGCAGCAAAGTGAATTTCTGTCCGGCTGCAATTTAAAAATGAAAACCT AGGTAGATTGACTCTTCTGCATGTTTCTCATAGAGCAGAAAAGTCTAATCATTAGCCACTTAGTGAT GTAAGCAAGAAGCATAGGAGATAAAACCCCACTGAGATGCCTCTCATGCCTCAGCTGGGACCCACCGT GTAGACACACGACATGCAAGAGTTGCAGCGGCTGCTCCAACCTCACTGCTCACCTCTTCTGTGAGCAGG AAAAGAACCCTACTGACATGCATGGTTAACTTCTCATCAGAACTCTGCCCTTCTTCTGTTCTTTTG TGCTTTCAAATAACTAACACGAACCTCCAGAAAATTAACATTTGAACTTAGCTGTAACTTAAACTGAC CTTTCCCGTACTAACGTTTGGTTTCCCGTGTGGCATGTTTTCTGAGCGTTCCTACTTTAAAGCATGG AACATGCAGGTGATTTGGGAAGGTAGAAAAGACCTGAGAAAACGAGCCTGTTTCAGAGGAACATCGTCA CAACGAATACTTCTGGAAGCTTAACAAAACCTGCTGTCCTTTTTATTGTTTTTAATTAATATTT TTGTTTTAATTGATAGCAAAAATAGTTTATGGGTTTGGAACTTGCATGAAAATATTTAGCCCCCTCAG ATGTTCTGCAGTGCTGAAATTCATCCTACGGAAGTAACCGCAAACTCTAGAGGGGGAGTTGAGCAGG CGCCAGGGCTGTCATCAACATGGATATGACATTTACAAACAGTGACTAGTTGAATCCCTTGTAAACGTAG TAGTTGTCTGCTCTTTGTCCATGTGTTAATGAGGACTGCAAAGTCCCTTCTGTTGTGATTCTTAGGACT TTTCTCAAGAGGAAATCTGGATTTCCACCTACCGCTTACCTGAAATGCAGGATCACCTACTTACTGTA TTCTACATTATTATGACATAGTATAATGAGACAATATCAAAGTAAACATGTAATGACAATACATAC TAACATTCTTGTAGGAGTGGTTAGAGAAGCTGATGCCTCATTTCTACATTCTGTCATTAGCTATTATCA TCTAACGTTTCAGTGATCCTTACAGAAAATAAGCAGCATATGAATAACCTGCAAAAAAAAAAAAAAAAA AA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG</pre>
<b>Restriction Sites:</b>	Sgfl-MluI
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
<b>Components:</b>	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
<b>RefSeq:</b>	<a href="#">NM_001159704.1</a>
<b>Summary:</b>	This gene encodes a member of the four-and-a-half-LIM-only protein family. Family members contain two highly conserved, tandemly arranged, zinc finger domains with four highly conserved cysteines binding a zinc atom in each zinc finger. Expression of these family members occurs in a cell- and tissue-specific mode and these proteins are involved in many cellular processes. Mutations in this gene have been found in patients with Emery-Dreifuss muscular dystrophy. Multiple alternately spliced transcript variants which encode different protein isoforms have been described.[provided by RefSeq, Nov 2009]

Locus ID: 2273

MW: 52.6