

Product datasheet for **SC214177**

CASKIN1 (NM_020764) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CASKIN1 (NM_020764) Human 3' UTR Clone
Symbol:	CASKIN1
Synonyms:	ANKS5A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020764
Insert Size:	1400 bp



[View online »](#)

Insert Sequence:	>SC214177 3'UTR clone of NM_020764 The sequence shown below is from the reference sequence of NM_020764. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GCCGACCAGCTGGATGCCATGCTGGAGTGAACGCCGCTGGCCGGGCCCTCCCGCGCCGCCCGGGCCCT CCCCGCACACTGACCTATACCTCAGGATGGGCGCGTCTGGGCGCGGCGGAGCGGCCGGGAGGGCCCT GCAGAAGCACAACCTCCGGCCCTGGACCCGGGCGGGCGCCACCGGGGACCGCTGGCCGGGGCTC CAAGGGCTCTAGGCAGACCCTGCGCCCGCGGGTCTCCCACTGCTGCTCCCGTGCAATACCTGCTG GGCCTCTGCCCCGCGGGAGGGGAGGGCGCCCGGGACCAGGGATGGGCGAGCGCACAGGCCCGGGCC CAGCACAGAACTGCCATCGGGACGCCCGGCCAGCGCTCGGGCGCAGCACAAGGACAGGGGACCAAG GTCAGGGCCCCCTGCCACCCGACCCCCAGAATATAAGCTATCAAGAGTATTAATTTATTGGGAATG AGCTGAGCGGATTTCCCAGAGAAAACAAAAGGTATAACTTTAACAAATATATATTTAAAGAAAGAAA TATTATTGATTCTATAGAAAACCATTTACCAACTGAAAGGACACGAAGCTAGCTCCGGACAAAAGCTG TTAGAGGCCCAGGCTGTCTGCCTGCGGTCTTCTCCGGCCAGAGTGGAGGCCCAAGGCTGCTAGCACA GGTGGAAGGCTGGGGCCCGGGGCGGGTGTCTGCGGGGGCTCTCCGTGCTTCTCTCCTAAGTTCTC TGCGCCCCAAGAACTCGAGGTTGTGCCTTGCTCTCGGCGGCCGTGTCTCTCTCTGGCGTGGCTGCGTC CCAGGGCGCGGGTGGGTTTGGTGAATGTGGTGTGTGTCAGCGTGGCAGCCACGTGGAAGCAGGGCATGTA GCAGGCCAGGCCACGGGCTTGACGCTCCTCTGAGGGTGTGCCTGGCCCCCTCTCCAGGGGCTGTGCAGG CAGCCCATCCACCCAGTGGAGTTCCCACGGACGCTTCTGCTCATGGTCCAAGGCTGTGTGTGTGGTC ATTCTGCCCCCTCCCTTTCTCTCTCTCTGCCCCCTCAAGAAGTACTGACTGGTCTCTCTCTATCGT GTTTGCTGATCCACAGTGTGGGCACAACTTTGACGTTTCTTAAAAAATTTAAAAAATGCCGAAA ACCAGCGATTGTGTGTGGGCGGGTGGAGGCGGTGGCGGAGGTGGCGTCCCGCGTGCAGGGCCCT TCCTGCTGTCTCACCCGTGTTCCCTGCGAGTGGACGTGCCAGCTGTCCCCCAGGCCTCTCTCAG AGGGCGGCACAGGCTGGGCTCCAGCAGGTGGGAGGGAGGTCTGGCAAGGCGTGCCTGCGGTGTC AATAACAGCAGGAAACAGA ACGCGTAAGCGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_020764.4
Summary:	May link the scaffolding protein CASK to downstream intracellular effectors.[UniProtKB/Swiss-Prot Function]
Locus ID:	57524
MW:	49.2