

Product datasheet for **SC214179**

ZNFX1 (NM_021035) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNFX1 (NM_021035) Human 3' UTR Clone
Symbol:	ZNFX1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_021035
Insert Size:	1397 bp



[View online »](#)

Insert Sequence: >SC214179 3'UTR clone of NM_021035
The sequence shown below is from the reference sequence of NM_021035. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AACTTTGAGGAGATCCAGGGGATGATAGGAAGATGGTACACCACTGCCTTTTGCCTCGCCACTGAA
TGACTGGGGCCAGCTCCCTAATGAAGGAACTGAAGTTTGTATTTTATTATCATCCTTTTATAGGCTGGGC
GCAGTGGCTTACGCCTGTAATCCAGCACTTTGGGAGGCCGAGGCAGGCGGATCACGAGGTCAGGAGTT
CGAGACCAGCCTGACCAACATGGCGAAACCCCGTCTCTACTAAAAATACAAAAATTAGCTGGGCGTTAT
GGCGGGCGCCTGTAATCCAGCTACTTGGGAGGCTGAGGCAGAAGAATCGCTTAAACCCAGGAGGCGGA
GGTTGCAGTGAGCTGAGATCATGCCATTGCACTCCAGTCTGGGCGACAGGAGCAAGACTCTGTCTCAAA
AAAAAAAAAATCATTCTTTTAGTCTTAGCACCTACTTAAGGATCCACTTTTAGGGCTCACCCACATTT
GTTTCTAGATTTACCCCTGCGCTAGAGTAAGCACTTTATCTCCAGAACTGAGAGCAAAGTTAACAAATC
TACCCCTTCTCTCTGCAAATTAGTGGACAGACTCCCTGGAACATGTTTGGGGCTTCCACCTAGGGCC
ACCTAGTGGTATCTCTGGGTCTTTACTTGGTCAGATGTTTATTCTACATTGTTCCCAAGAACAGAGTA
TGAGCTCATTGATGCAGACCGATTCTAATTGCCAGGCCCTAATTGCAGACTAACTCTCATAATAAACA
GAGGCCATAGTTGTTTATGAAGTGCTTATCCCTTAAAGGAGCACAAGAACCCTCCCTGCCCTCCTTG
GGCACCCTGCCTCCAGGAGATGGAGGCACGTGATAAGACAAAAGACTGCACCAACTCACCTGACACAG
TTACATAGTCACTGAGAGTGGGAAGATGGGACAGCCACATGCTGCATAAGATGGGCCTTATGCAGCA
GGCCCAGGTCGTCATTAAAGGAGTGACCCCTTTCCTGTAACCTGCACTTTGGGATGGTAGAAGTTTCTTT
ACCTGCTGACAGGTTTGGTGGCACTGCTGTTACCCCTGGGCCCTGAATGGAGCTAAAATCACATTTGG
TACCAGCAGCACCTATCCCAAGTGATCCTTCATCCCAACACTCCCTCTTGGAGCTGTTCCCTGGGTA
GAGCTAGCATGCCAGCAGCTTCTGCAGGCTCCAACCCAGGCCAGAAGCCAGACCCAGGCCTGCTGCCT
GCATCTGCATTCCCTCCTTCCAGTGTTCTTAGAACAGACATTTAGGTATCTCAGGTCTTTCTAAGTG
TCCCTTTCCTATGTATGCATTTCTTTTTTGTCTTTACTATGCACTTTAGCTTATAAAGCCAATTA
AACATGATTGAGAAAA
ACGCGTAAGCGCCGCGGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_021035.3](#)

Locus ID: 57169

MW: 51.9