

Product datasheet for **SC214181**

C6orf64 (SAYSD1) (NM_018322) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C6orf64 (SAYSD1) (NM_018322) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SAYSD1
Synonyms:	C6orf64
ACCN:	NM_018322
Insert Size:	1399 bp



[View online »](#)

Insert Sequence: >SC214181 3'UTR clone of NM_018322
 The sequence shown below is from the reference sequence of NM_018322. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTACAGTTGAGACCCCTGGCAGGGAGATAGGACCCAGCTGTGCTGTCATGCAGCTAACCTCTGATGTGG
TCTTCCTCACCATTGGCTATGGATTTGATTTTCAGGTGTATAGGACTAAGGGCAGCTTGCGGGTTAGCTC
TGTGACTGCATAGTTTTTCTACCTTCTTTCCCTGATCTTTTGTGCCATTTGATCTTTGATAGTTTTGG
TGAAACTCTCTAAAATACATTCACTGTGGGTCGACGCAATTTATAAAAAATTATGTAACAAGAAGGGA
GACCTGTTTGTTCATTTCTCATCTGTTTGGGAGATGATTTTAGAGCACTAGAAAGGCACTGGGGAGAT
TCTCAGCTTAAAACATCCAGCAGTTTGAAGTATGATTAGGTACATCAGGGCTGCATTGTCAATGTTCTC
TTTAAAGTCTTTTAAACATTTATAGCAATTTTTTTTTTCCCGGAGAGTTAGGTTGCAAGTTTTGGGTTTC
TTGTTTGTGTTTTGTTTGTCTCCTGCTTAATTCTTTAATTTTTCAGTCATTACTGGTATTGAAAAATAA
AATATCTTTAAAACATCTCCTCTCAGAAATAGGTCCTCTTCATTGCCATCACCATCTCCCACTCTC
CTATTATTTTGGCACTACTCAGTAAAGGAAGGTAGGAAGAGACAAACGCCTAAGTGCAGGTGTGGGGAG
GGATTTACAAGTGGTTATTAACGGCCAGTTCAGCAAGAAGTGTGAGTGTGTACAAAGGGGAGGGCTG
GAAGTGTAACTCCAGACCCGTTGGCTGCTTGAGTTGTTTCTTATATTCTAAAGCAGCAGTCCCTAACCC
TTTTTGGCACCAGGACCCAGTTTTGTGGAACACAGTTTTTCCATGGACGGGGTGGTGGTGGAGGATGAA
ACTTCCACCTCAGATCATCAGGCATTAGAGTCTCATAAGGAGCACGCAACCTAGATCCCTCGCATGCGC
AGTTCACAATACGGTTCTAAGGGCTTTAGAGTAAGCAGCTTTTTCACCTGTGGCCTCTGGTGAGAAAT
TCTGTAATTTGTGATAATCAGGCTGGATTTAATGCTGCTTTTCCAGTACAATGTTAGAGTTTGGGTTTC
ATTAATAATTAGGCAAACTCCCATTTGGGTTAGGGCTTCTCTCATTCCATTTTGTGGCTAACCTTACTGTG
TTTCAGCCCTTGCTGAAAATTCTCTGATATGTGTTGCCCTTCTCACAGCCCTTTGGCCATTGGGAGT
TTGGCTGTCCCTCAGAGCCATCCGGTCAAGCAGATGGTCTGTTCTATCTCACAGAAAAGTCTTTTCTTC
CATGAGTTCTGTCTGAATGAACATGTA AAAAGTATGGGAAACAGATGAATCCCTATTAACATGAAGT
TTTGATTGTATTTAAGATA
ACGCGT AAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_018322.3](#)

Locus ID: 55776

MW: 52