

Product datasheet for **SC214192**

SCYL1BP1 (GORAB) (NM_152281) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SCYL1BP1 (GORAB) (NM_152281) Human 3' UTR Clone
Symbol:	SCYL1BP1
Synonyms:	GO; NTKLBP1; SCYL1BP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152281
Insert Size:	1402 bp



[View online »](#)

Insert Sequence: >SC214192 3'UTR clone of NM_152281
The sequence shown below is from the reference sequence of NM_152281. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AATGACATTTTCAGCTGCTTTGGCCACATGAAGTTCTGGTATTCTTTTGGAGCTAATATGGTATTGAGTAA
AGTATACTTTTTGCAGTAGATCATGCCCTGACCTCCAATAAAAAACCTCTTTAAAAACAATGCTGATTTTT
GAATTCATGATTAGTTTTAGTAAATTAATAGGTTTGGCCATTCTAAAATCCAAAATTATGATAACTGTA
TTTCTAGGGGTATGTTTCACCTTGATTTTGGCCCGTTCTTTTCAGTGTTCCGTTTACCCTTTTACTGA
AAGTAAGTGACTTAAAGGATGGAAGAAAACTACATGGTTGGAGTGTTTTAGATCAGATAAATATAG
AAAAGTAAATGGAAATGATACTCAATACCTAGTTAACTCATTTCACTTCTCTGGGAGACGTATTA
AGACTGGAAGTCTATTTTATTTAATCAGTATCAGTAAAAAGACAAAAATTACAGTGTGAACAAGAGAA
AAACAGGAATAGAAGAAAAAGTTGCAGGTTGAATTATCTATCTGGATATTACTAGAGTTGTAAACAC
AACTTGGAAACTGGAATTTATGTGAAGAACCAAACTTAAACCAGCATCACTTTTGTCTTCAATTTG
CCTTCCAGTGATTTTAGTTTACACCTGCGTTTGTATTGAAACACACACACTTATATATAAACACA
TGTAATTTTACTATTTTATTGCTTTTCTTTAATCGATGAATTGATTAATTTAGACAATTAACATT
TCTAAAGTGAGACTGAAAATAATATAGAAAATCTCATTGCTACTTGTGATTATCAAAAAAGTATGTAT
TTTCTTTCTAAAGTTATCTATAAATATGTTTTTGATATGTTGGACATGAATTACCCTGTCATGAAGCGT
TTAAATTTGTTAAATGCTGATTTTCTTATTAGTTTGTTAATTGCCTGTTACTAGTACTTGACCAATACA
TTGTGCTATATGTAATGTGTAATAAATTGTAATTTTCTCACAAAGTCTGTGAATATCACATTATG
ATTTATTTATCAGAAATCCAAAAAGTAAAAATATTACCTTTGAGGACTTTTTTTTTTAAAAAAGAA
TTATTATCAGAGAACATAAGCACAAGTATAAGATTGTATCCATAAAAGTTTCTGAAAAACTTAGGGACA
TCCTGTTTTTAAATGGGTTTGTGTTGAGGGTATTATCTTCTGGTAATATAGATGAGCTTTATATTTTAT
ATTATGAATTTGAGTTGAAAGTATAACATTTTACTACTGACTTCTTAAACAGAAGCATTTCTATAGA
TAGTTGTGCTTTTATTTATGAAATACATAAAAGCCAACAGTTTCATTGTCTAACACACTTCTTTTCTA
AGGAATAAACAAGTCTGATGTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_152281.3](#)

Summary: This gene encodes a member of the golgin family, a group of coiled-coil proteins localized to the Golgi. The encoded protein may function in the secretory pathway. The encoded protein, which also localizes to the cytoplasm, was identified by interactions with the N-terminal kinase-like protein, and thus it may function in mitosis. Mutations in this gene have been associated with geroderma osteodysplastica. Alternatively spliced transcript variants have been described. [provided by RefSeq, Mar 2009]

Locus ID: 92344

MW: 54.7