

Product datasheet for **SC214199**

UQCC (UQCC1) (NM_199487) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	UQCC (UQCC1) (NM_199487) Human 3' UTR Clone
Symbol:	UQCC
Synonyms:	BFZB; C20orf44; CBP3; UQCC
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_199487
Insert Size:	1387 bp



[View online »](#)

Insert Sequence: >SC214199 3'UTR clone of NM_199487
 The sequence shown below is from the reference sequence of NM_199487. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCTCCGACTTACAACGACGAGGGACTTTGATGGGCTGGGCCCTCCGCACGGCCCGCCAGCTGGCTTCGA
GGAACCTCCAGGAGAGAAGTGCCTGTTGGTCCAGGACCTGCAGAAAGTGGCCTGAACTGACCTCTGAA
CAGCATCTGTCAAATACCTGGCCCCATTGTGTTGAGTTTCTCTTAGTGTGCCAGGAGTCTGATCTG
CTGGGGTACAGGGCTGGGAGAACCCTAGCTCTCCGGGGTGTCTCTCCCTTAGGGGAAGCCCCGAGT
GAGAGTCCCCAGCACACTCCCCAACCCCTCCAGCAACTACATGTGACTGATAGCTTTTCCAAAG
GCCAAGGAAGGGATGGTGTAGGTTCAAAGGGAAACCCCCAGGGCTGTGTGGCTAGGAGCAGATT
GTAATGCTGCCGAGTCCGGTCGGTGACCACGCTTGTCCCTCGGCTTTCAGCCATGGGGTTGAGTTGGC
CATTAAAAGAAACAGAGACTTCTCTGCCATGGCCCTTTTATCCAGGGACTTAGAAACTTGCCTG
AGATGGTGGACGAGTAATGAGGGCACCCGCGCAGCTCAGTTAGAGACGGAGAAAGGGAAGAGGCTGGGA
TGGTCTCTGCTGCTCTTGCCCTAGTTCATGGAGATGTGTCTCTGTTCAGGCCAAGATACAGCCAGCCA
GGCCTGTGCTGCTGGGACCCAGGAGGCCTCTGATGACCAAGGGCTTTCACATCCTAAGTCAATTTGGAAGG
AGGCCTTGAGAACAAAGTCACTTTGTCCTCCAGTGAAGTGAAGGAAACATGCTGTCTCTGTCT
TGGCCTCCCTTTTATGAGATACTGGGGAGAAGAGAACATTCTCCTGGCTTAGTTGTAGCAGACCCAG
ACCTGTGCCAGCTTTGGTCCCCCTTCCAACTTCTGAAGCACGTGCTGCAGAGCCACCTTGGTCTGAG
CACCTGAGGACCAGCCCTCCTCCCTCAGTGGGGTTCATCTCTGGGGGATTTTCTTAAAGTGAAGAAA
GGGGGTGGGGAACCATATTGCCCTCCCTCCCCATCAAACCTTCTTAACTTACTATAAAATGA
GTATATAAAGAACTCTATATGGGTGAGGTATATCCCACTTCTGTGAAAACATTACAAATCAAACCGC
TTCTCTCAGTTTTATTAAGATGCTTTTGTGCGAGCGGAGCTCTAGAGTGAAGCCTCTGTGTGTGTG
GAGATAATAACACCTTGAACCTATTACAGCTGGGCACTATTTACATAAACCCAGAGCTGAGCCAGGCAG
GAATTTGCTGATTAATTTATTTAATGGAGTGAAGTATACCATGCACCAAAAATAAACTTTACTGTGTG
TACCTAA
AGCGGACCGACTTACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCC
CAACCTGCCATCAGAGATTTGATTCCACCGCCGC
  
```

Restriction Sites: Sgfl-RsrII

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_199487.3](#)

Summary: This gene encodes a transmembrane protein that is structurally similar to the mouse basic fibroblast growth factor repressed ZIC-binding protein. In mouse this protein may be involved in fibroblast growth factor regulated growth control. In humans, polymorphisms in this gene are associated with variation in human height and osteoarthritis. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]

Locus ID: 55245

MW: 51.1