

Product datasheet for **SC214228**

TENT2 (NM_001114393) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TENT2 (NM_001114393) Human 3' UTR Clone
Symbol:	TENT2
Synonyms:	APD4; GLD2; PAPD4; TUT2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001114393
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC214228 3'UTR clone of NM_001114393
The sequence shown below is from the reference sequence of NM_001114393. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

CGACAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGCCGGAAGATCGCCGT
TAACAATTGGCAGAGCTCAGAATTCAA**CGCATCGC**
CCTGTAAGAGCTGCTGCTCCTGAAAAGAT**TA**ACTGGCCTCTATTTCTTAATAAATCTTCCAAGAAATAAA
GAACAATAGTTTCATCATAATACATTATGTTTACCTCCATCATAGTTGCTTTTTTCATAGTTCTTGTTT
TCATGTTTTATTTTTAAAAAGACATATAAAGATTGCATATTTTATTATGACATTTCTAATAAAACCCCT
TTGATTTAAAAATGTATTTTTAAAAATGTTTACATTGATGATGAGTAATATTGCATGTGTTTTCAGGTG
ATCAGATAAAATATATTTTGGGAAATTAACGAACAAATAAAAAAGTTTTTGATATAACTTCAATTAAT
GTACCACATGCTAATCCTGAAGAGATGTGTAGAATTTTGGAAAGGCTAATTCATACATGCTTAAAGT
ATAACCTACGTTAGATAGTTTGTGTCAAAGTCTGTTCTAGTAAAGTGAAGATTCAAGTCAGTTGTTCA
GTTACTTGAAGCAAAACGAAATCTTTCATTTCAAGTCAAATCACTGCAGTCATGAAATACTGAACAATTG
CCTTAAGTCTTTGCTTGACTCACTGGGATAGACTGAGGCTTTGGGTGTGCTGTATTAGCATTTTCATTA
GTACTTCACATGCTTTTTGATGTACTCTTGAGATTGCTTTAAATTTGTATTGAAACAACAATACATTTT
GCAGTGTAGTAATGGGAGCACAACCTCTAACACAGTTAGTGAATCGTTTTAAAGAATCAGTTTCAGTGT
GACATTTTTGAAAAGATTGTTTCTGTGCTTTACGATAGCTTAGTGCAATGTGCACCTTCTGTTTTACTT
GCCATTTTCTGTCTGTGTTTTCTGTGTGACATGAAGCAACAGAACTGAGATCAAAGTTAAGATTATAT
CCTGTTTGTAGTATCAGATATTTTTCTGTGTACAATTATAGGATTGTAATCTAAACTGGAATTTTTAGG
CAGTAAGTCACCACAAAAATGTTTTAGATAAGACACAATAAAATTATTATAAATAAAAGCTTAATGTTTG
TAAAAAATCTCTTTTTTAGTATTTCTTTTTTCACATGAAAGAAGTGGTGGCTGCTAAAAAAAAGCTAC
AGTGTTTATTAAGGGTCTTTTTGATTTATGTAATATTTGTAATTTGGTCAGTGCCTGTAAATTTAAAT
ATAAAAAGTAACCTTGAAAACAGTTTTAACTTTTTCAAAGAAGTATGTCCAACATTTTTAGACCTGC
TGTAGTACAGTTTTGTACCTCTAACGTATTTTTTTTTTGACAGCCAAATGCTAAACTTTTTGCTTTTCT
TTGACTTGTAAGAAGTGCACATTTTCATTTCTTCCCTAAGTTCAAATTTTTGTATGATGTCAAATGCA
ATAAAATTTATATATGGACATTGTTGAGGTGCCTTTTTTTTTCTGTGAAAATCTTTGCTTTGAGATGAGA
CCGAGGCTGTAGTTATCATCTAGAAGTGAAGTAAATCTACTTATTAGCCAGCTAACGACATCAATATCT
TATGCTCTGGCATTATGAGGAATGAAAGTTTACACATAGTGAAATTTTTTAGATTCTGATGTGTAACA
CTTTTTATGCATTTGAGCTTTTCATATTTTGGGGGGGAATTTTTCTTAAATTTACACTGTGCATCTCT
TAATCAGAATGTACTGAGTAGAGTTAATTTCTAGTGATGACTGGGTAAATCATTATAAGCATTAGTCAT
AGTACTATTTCAGTAATACTGTGCAACTTTACGGGATGGGTAAAGTATGTACCCTCTGCCTCTACAGTATT
ATGACATTTTCTTAATAAAAACCCAGAATTTACCATGGTGTGGGTTCTCTACTGTCACTGTTTCTCAG
TTTTTCAGTGGCCCCACACTTACTACTGCTTTTTCTTTGTTTTACCACGAAAAACAGATCAAGAAGTCACT
AAACTTAGACCCAATTTCTATATCAGTAGACTGAGGCTTGAATTTATGCACACATAACTGGCCAAATAT
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTGATTCCACCGCCGCTTCTATGAAAGG

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001114393.3](#)

Summary:	Cytoplasmic poly(A) RNA polymerase that adds successive AMP monomers to the 3'-end of specific RNAs, forming a poly(A) tail. In contrast to the canonical nuclear poly(A) RNA polymerase, it only adds poly(A) to selected cytoplasmic mRNAs (PubMed:15070731). Does not play a role in replication-dependent histone mRNA degradation (PubMed:18172165). Adds a single nucleotide to the 3' end of specific miRNAs, monoadenylation stabilizes and prolongs the activity of some but not all miRNAs (PubMed:23200856).[UniProtKB/Swiss-Prot Function]
Locus ID:	167153
MW:	77.3