

Product datasheet for **SC214231**

USP28 (NM_020886) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	USP28 (NM_020886) Human 3' UTR Clone
Symbol:	USP28
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_020886
Insert Size:	1395 bp



[View online »](#)

Insert Sequence:	>SC214231 3'UTR clone of NM_020886 The sequence shown below is from the reference sequence of NM_020886. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CAGGGAGTTTCAACTGTGACAGTGAAATAAGTCTCCACATGTTCAAGGCCATTCTGGTTCCTGGCTGC CTGCCTCTTGACAGAAGTTCGTTGTGATAGTGCTCACCTTGGGAAAAGGATTAGGTGGGCACATAAGA TTCCGATCAGACCCCAACCATGCTGCATGTGTAAGAAGGATTGAAAATAAAATTGCACTTTTTAGGTA CAAAATCATAAAAGCTGTTTCACTAGAAAAGGCAGAAAGCAGTGTATTAAGGTGTTGAATTACGCCAGA AGACCTGAAATGCCTTGACCTACAACAATGCTTAGGCTTTTCTAAGCCTCTTGCCACTTTTAAATTA TCCTTCAGGCATAAATATTTTTGACAGCAGAATAGAAGAATGATTGATGAGAACCTGAACCAGATGAAC AGCTACTAGTTATTTTATCAAAACAGATGACATTTAAAAATTCTTAACAAGAGATTAGAAATATA AACCTTGCTGGCTCTTGCCAGGAGATAACAAAATGGGTTGCTGATGAACTGCACCTTTTACATGTGG GTAGAATATAAGCTCACATGGCAGTGAGATGTTGAAAAGTCAAAAGAGACCTGTCTCTCTCTTTCTTT TCTATCTTTAAACCAGAAAACCTCATACTCAGTCTCAGTGAAAGAAAGTAAAGTATTAAGGACTTTAG ACAGAAGAGCATTGTGTAACCTTGACTGAAGATCATCCATTAATAGTTATTAGGCATTTAGGTAAATTT TCTAATACCTAAAAATTGTCAAAAACAGTCAATAGGGCTACTGCTGGCCCAAAGACCATTTAGGTCCAC CTCCTCTTTTTTGTCTTTTTTTTTTCTGTGACAGTTTCACTGTGTCGCCCAGGCTGGCGTTCAAGTGG TGCAATCTCAGCTCACTGCAACTCTGTCTCCTGGGCTCAAGTGATTCTCGTGCCTCAGCCTCCCGAAT AGCTGGAATTACGGGCATGCACCACACACCTGGCTAATTTTTGTATTTTAATAGAGATGGGGTTTCA CCATATTGGCCAGGCTGATCTCTAACTCCTGGCCTCAAGTGATCTATCTGCCTCCCTCAGCCTCCCAAA GTCTGGGATTGCAGACAAGTCATCGTACCGGCCCTTCTTTTTTGGCCTTAAAGTAAGGGATGTGGGTT TGTACAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAGAAAAAAG AGAGAAAAAATACTTCCATTGATCATTTTTTAAAGGCTTCTGATTGGATATTGTGTTTTAACCAATTT TAAAGATTAATGGAATCATGAAAGGGAAAAAATTGATACAACATATGCAGATTTTATAAATGTGCAATA AAAGTATTTGTTTTA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_020886.4
Summary:	The protein encoded by this gene is a deubiquitinase involved in the DNA damage pathway and DNA damage-induced apoptosis. Overexpression of this gene is seen in several cancers. [provided by RefSeq, Oct 2016]
Locus ID:	57646
MW:	53.5