

Product datasheet for **SC214243**

TSPAN17 (NM_012171) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	TSPAN17 (NM_012171) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	TSPAN17
Synonyms:	FBX23; FBXO23; TM4SF17
ACCN:	NM_012171
Insert Size:	1383 bp



[View online »](#)

Insert Sequence:	<p>>SC214243 3'UTR clone of NM_012171 The sequence shown below is from the reference sequence of NM_012171. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC TACCCTGAGCCAACCTTTAAAAATTGGTAGATTTACATAAAAAGTCCAGATCCACAGCTTCTCTTGAAG AATGACCACCTGGCTACGCCGGCTCTTCGGTGGCAAACTACCTGGGACACTGCCTCCCCAGTCACCAA GGGCCCCAGCTGGCCCGTTCTACTCACCTAAGTCCGCTGACCTTGTACACTAGGAGCTGGCCTCCC ACCTCTGCAGGGTTATTCCTGCACCTCGAGGCCGCTGCGGGCCAATCTGGAGTGAACACGGGGACCT GAAGGATGGAGAGGCTGGACCCCGCTTTGAAGAGGGTGCAGCCTGGGAAGGGCGGCTTGTGGGGACT GCGGTGGGAGTAGAGTCCCAGGAGAGGGTCTGAGGGTGGGATGGGGTTCAGGACAATTTGCAAAAG AAGTAGCTGGAAGCCATGGGACTGGCGGGAGCCTGTTGGGGATCTGGATGGTTGACTCCTAGGAGTC AAGTTCAGCATCTTCACCGTGGCTGCAGAGCTGCCTGATGGGCACTAGAGGGCATGCCAGCCCCACT CCCTGGGTCTGGCTTCTCCCGCAACCTCACTCTAGTAGAGCCTGTGCCTACTAGCGCTCTGGGG TTCGGAGAGTTTGGGAATTTCTCAGAGCCAACCTGGCTCAGGCTTGGGAAGGCTGGCTGCTGCCCTCAGC TCCGCCTCATCAGCTATGTGAAGGGTGTGTGGAGTGATCCTGCCGCCCTCCCTGGGCTCGTCCA GAGATCTCAAACCTCCGATGCCCTGGGGCCACGTATGTTGTATAAATGGATGAAACAGGCCCTTGAGTT GGGAGCCTGCTTCACTTTGACTTTCCACTGTTGCTGGAGACAAAGACATCGTGATGAGAGAAAGTTG CACAATCTAGTCGGTAACAGCCACTTTCTTGAGACCAAGAGAGTGCGGTGGGGATGGGGGGGAGAGCA CGGGTCCCCTGTCGACAGTGGCCGCTGCCATATTCAGGTGTAGCTAATTGCTCTGGTGTGGGAATGCAG GCCTAATGACAGAAATCTGGAGAAGCCAGAAATACAGATTTGTATGTGAGATGTCTGATTTTTTAAGT TGTTGGCAGAAATTAATTCAGAAATCAAATCTGCAGGCCAAACAAGGTGCAGGCCAGCTTTGGCCCC ATGCCCTGTAGTCCCTCTGGGACAGTACCCTGGGGTCTGGCTGCTCTGTCATTGAGGGATGCTG GGCACTGCTGCCGGTGGCCAGGGTATGGGGCATGTGCCAGCAATGTGGCTCCTTGCCCGCTGGCC AGTGTCTGGGCCCTGACAGGGCTGGCTGTGAGTGGTTGTACATGCTACAATAAATGCAGCTGGCA GCA ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_012171.3
Summary:	This gene encodes a member of the transmembrane 4 superfamily. It is characterized by four tetraspanin transmembrane segments. The function of this gene has not yet been determined. [provided by RefSeq, Mar 2014]
Locus ID:	26262
MW:	48.8