

Product datasheet for **SC214270**

SR1 (SRI) (NM_198901) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SR1 (SRI) (NM_198901) Human 3' UTR Clone
Symbol:	SR1
Synonyms:	CP-22; CP22; SCN; V19
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_198901
Insert Size:	1393 bp



[View online »](#)

Insert Sequence: >SC214270 3'UTR clone of NM_198901
The sequence shown below is from the reference sequence of NM_198901. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GATTTTCATTCAATGTGTCATGAGTGTTAAATCAAGAGGAAGCTGCATGAATGTAATCAACATTC AAC
TGGAGCTCTCCTTTGCTTGCCTCTTTGCCTTCGGTAATATGTATAAACTTACATCACGACTTTCTCTT
AACAGCTGTTGTAAGTTTATTACTTTATGTACAACCTGAAGTTTTGTTTTAGTTTTGATAATAAATTCT
TTGGAACTTTAATAAGATCTAGTCTGTTACACCATTTAGAACTTCTGAGCCATTATCAGTCATGCCT
TATTTTCTTGCTAAAACCTATGTAATTTAAGTATGCAAAATGTTTAAGTCACATTATTTATTTTCA
TTGTGAGACACTAAAACCTGTTAATCAGACTACAGCTGTTATCTTCTCTCCTACAAAAGAATACTCCA
CACATAAAAACCTTAGGTAATGACATAGACGCACTTGGGTGAAATAAAAACAACAAAAAGGTAATCCAG
TAATCCACGTCAGGATTCACCTTAGAAGTTTAGCACACGCCCTTCAAACCTGTTGAATAATTTGATTG
GCAAACTATCTGTCCCAAGTCCCTTTTTGTATCTATTTAAACCTTTGTTAACTCTCCTTAAAAAT
CTTGATACATTATAAGCTTAACTATATAAAAAGAAAATTGATAGAAAAGACTAAGGGAATATATGGAAA
TATTAACAATGTTTATTTGGATGAATGGAATTCTAGATGACTTTAATTTCTTTATATCAGTATTTTT
CAGATTTTCTATATGCTTTCTTTTAAAATCAGATTTTAGTTTTTAAAATACTGATACATCCATTTTGA
TATACCATGTATTCACACTTTAAAAATACACATAGCAGAGTATCGGAAAGTATGTACAAAAACATTGA
TCTGGCTACCACTGGGAGTCAGGATCATGAGTTCATCTTCTTCTTCTTACTGCTTTTCCAAACATTCT
CCAGTAAGCAGGTAACATTTAATGGAAGGAATTTTTTAAAAATTTTAAGCTGTACACTTTCAAG
TAAGATCTGAATTCTAATGCTGGCTTGTGCCTTACTATGTGGTTTGGTCATTATAGATAATGCCAG
TTTCAGTTTTCCCATCTTTGAAAAGGAGATAAATATGTTATCTCTTGGGAGTAGTCTCTGAATGAAGTG
TGTTGATACAAGAAGCACTGACAATGTTTCTGTTTCTTAAATTTAAAACCTGGCCTGGTTTGCCTTTTTT
ATCAAGAGAGCTTAACAGATAAAAAATGAAATTAGTCTATTTTCTACTTGCCAGCAGAGTATCTGTCTT
ATTTTAGGATGCAGTGTGAAACTTACCATTCAACTGAACAACCTAGTTGTCAACTAATGAAAAATAAACA
TTTTCTCTAAGCA
ACGCGTAAGCGGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_198901.2](#)

Summary: This gene encodes a calcium-binding protein with multiple E-F hand domains that relocates from the cytoplasm to the sarcoplasmic reticulum in response to elevated calcium levels. In addition to regulating intracellular calcium homeostasis it also modulates excitation-contraction coupling in the heart. Alternative splicing results in multiple transcript variants encoding distinct proteins. Multiple pseudogenes exist for this gene. [provided by RefSeq, Mar 2012]

Locus ID: 6717

MW: 54.1