

Product datasheet for **SC214296**

ZNF547 (NM_173631) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF547 (NM_173631) Human 3' UTR Clone
Symbol:	ZNF547
Synonyms:	MIP-2A; SEDLP1; TRAPPC2.19; TRAPPC2B; TRAPPC2P1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_173631
Insert Size:	1410 bp



[View online »](#)

Insert Sequence: >SC214296 3'UTR clone of NM_173631
The sequence shown below is from the reference sequence of NM_173631. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTCAGACATCAGAAAGTCCACACTGGATAGGCCCTTATGAATGCAGTGGATATGGGAAAGCCTTCAGT
CACCAACATATTGTGGCTGGACAGCAGGAGTACACACTGGAGAAAGACTGAATGCCGTGAACGTGGGT
AATTATGTAGGTACAGCTCTCCAGTCGCTATGTATCAGAGAATTCACACTGCAGAAATGTGTGTTACGC
AAACTCGGACATTATTTTGGTTTGACTCTCATCTCATTAGACATTGGAGAGTTTACACTGAAGAAGAG
TCTTTTCAATAAAGTAGAAAGTGGTAAAGATTCAACATGCAAGATTGTACTTATTGGGCTTCAGAATAT
CCACACTAGTGAAAGTCTTCTGAGTACAGCAAATGTGTGACATTATTTTGTCTACTACTCCACACTACTT
AGACATCATGTAGTTCACACTGGAAAAAGGCCACGTATGTGCCTTGAATGTAGCCAAAATGACGAACAA
CACCCAGAAATCTGTGATTTAGCACTGAGAACTAGTATTATATGGTTTTTAAAAACAATGGTGAAGTA
CATGCCACATAAAATTTGCCATCTTAACCTATTGTAATGTCTTGTTAATACTTGAAGTACATTAACATT
GTTGAGCAAAGAATATCCTGAACTCTTTATCTTGTAAATGAACTCTATAACCACCATTAACAAAAACA
ACTCATTTCCACATTCTTCAGTCCCTGGCGACCACCATATTTTAAAGTCATTATGACTCTGACTATTCT
TGGTACTTCTCAGAAGAAGAATCTTAGAGTGTTTTTCATTATTTTTTCACTTAGGATAATATCCTCA
AAGTTCATTATTTTTAGCATGTGTCAGAAATTTAAGGCTGAGCAATATTTTATTGTTGCATTTAC
CACATTTCTTTATTTCTGACATCTATTCATGGACACTTGGGTACTTCTACCTCTGGCTATTGTGAAT
ATTGCTACTACAAGCATAAGGGCACAAATATCTCTTTGAGACCCTGCTTTCAATTCTTATTGGGTACCC
CAAAAAGTGGAAGTGTGATCAGTGGTACTTCTAACTTTAATCTTTGATGAAGTGAAGAACTTTTC
CATAAACTTGACCCGTTTTACATTCTTAACATTCCACAAAAGTTGTGATTTGTCCTCATTCTAATCCA
TAGAAGGCAACATTTCTAACAGTTTCTAGGTGATTTTGTGTTGCTGATGTAGGACTACATTTTGAG
AACTGTCTATCCATTTTATAACAATGCAGAGCTTCTCTGTAATACCCCTTGTCTATGTTCTAT
CCTGTGAGCCTTTTCAATTTGTTAATTTTCCAGTGTTTACAACCAGGAAATAGCTTGTCATCTTTGA
TACAAAACAGAAATCTTTGCAAAATACCAA
ACGCGTAAGCGCCGCGGCATCTAGATTCGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_173631.4](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 284306

MW: 54.7