

Product datasheet for **SC214297**

Methyltransferase like 6 (METTL6) (NM_152396) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: Methyltransferase like 6 (METTL6) (NM_152396) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: METTL6
ACCN: NM_152396
Insert Size: 1422 bp
Insert Sequence: >SC214297 3' UTR clone of NM_152396

The sequence shown below is from the reference sequence of NM_152396. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCA**ACGATCGC**

AGAACCCATCTCCTGTGGTCCTGGGCCTGGATCCTAAGTCC**TGA**CCTTTCATGAGGTTGGCATTATACC
TTCCCTTGAAGAGGAAAATTTAAATCACTCTTTCGTCTGGGTGCGGTGGCTCATGCCTGTAATCCCAGC
ACTTCGGGAGGCCGAGGTTGGGGATCATCTGAGGTCAAGGAGTTTGGAGCCAGCCTGGCCAACATGGTGA
GACCCCATCTCTACTAAAAATACAAAAAATTAGCTGGGCATAGTGGCGAGTGCCTGTAATCCCAGCTAC
TCAGGAGGCTGAGGCAAGAGAATCGCTTGAACCAAGAGGCGGAGACTGCAGTGAGCCAAGATCATGCCA
TTGCACTCCAGCCTGGGCGACAGAGTGAGACCCCTGTTTCATAAAATAAAATCACTCTTTTTACTTTA
TATATTTTTTGGAGACAGGATCTTGTGGTACGCAGGCTGGAGTGCAGTGGTGTGTGATCACAGCTCATTGT
AGCCTCGACCTCCCCTGGGCTCAAGCAGTCTTCTGCCTCAGCCTCCTGAGTAGCTAGGACTACAGGTGT
AAGCCACCATGCCTGGCCCTTAATTTTTTTAATGAGCATGTACAGTAGTCCCGCTTTATCCATGGTTTC
ACTTACCACAATTTCAATTACCTGTTTTCAACTGCAGTCAAGTACTTTTTCTTAAGATATTTTGGGGTTA
TAGGCTCATGCCTATAACCCAGCACTTTGAGAGGCCAAGGCCAGCAGATCCCTTGAGCCAGGAGTTCA
AGACCAGTGTGGCAACATGGTGAACCTCTGTCTACAGAAAATACAAAAATTAGCCAAATGTGGTGGC
ATGTGCCTGTAGTCCCAGCTGGGGTGGGAGTTAGCTTGGCCTGGAAGGTGGAGTTGCAGTAAGTCGA
GATCACGCTACTGCACTCCAGCCTGGGCGACAGAGCAAAAGTCTGTCTCAAAAAAACAATAAATTTG
AAAGATCACATTCACATAACTTTTATTACGGTATATTGTTAAAGTTGTTCTTTTTTTTAGAGTTTGGT
ACCATCTTGGTTTCAGGCAACCACTGGGGTCTTGAACATACCACCCGAGATACGGGGGGCTACTATG
TACTTTTTGTATTAATAAATAGGAAGTGTATGTGTGTGTGTTAGAGACATGTAATAAAATTTCAAGGCC
TGAAAAGTAGTCAAAGTTTTATAGAATAACATAATTTAGCCATATTGTTCACTGTAAAATTCACAAAGA
AAAGCCTAAGGACCTGTAGCACGTAAGAGATTCTTGTCTACTCATGTGACTGTAAGTGTCCCCCA
GGCAACATGTCAGACAACCATGGGGAAAAATGAGGCCAAGTAGACTAAAATAAGACACTCATATATTTT
CTAGCCTAAAGGAGTGTCTGT

ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCG



[View online »](#)

Restriction Sites:	Sgfi-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152396.2
Summary:	Probable methyltransferase.[UniProtKB/Swiss-Prot Function]
Locus ID:	131965