

Product datasheet for **SC214318**

ZNF418 (NM_133460) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF418 (NM_133460) Human 3' UTR Clone
Symbol:	ZNF418
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_133460
Insert Size:	1403 bp



[View online »](#)

Insert Sequence: >SC214318 3'UTR clone of NM_133460
The sequence shown below is from the reference sequence of NM_133460. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GTTACACAGAAAGAAGTCCTTACAAGTGAAAAGAAATTTGGAAATTTCTTAGCTAAACCTCTGTGCA
TCTTCTTGATCAGAGGGTCTTACTGGATCAGGACCTTATGAGTGTGACAAACGTGGGATATTCTTTAT
GCAGAAGTCTTGTATTACATACAGAAGAGCTCCCACTGCAGAAGGGCCTCTTGAGTGCATGAATG
TGAGAAAGCCTTCTGCCTTCTGTCATTGGATAACAGATTGTTCTCATAAGGAAACACTGTACACGTAC
AGGAAATATTATTCTTGTAAACATAACACTGGAGGAGATGCCTTATGACGGAGCCATCTGCCTAAAT
TGACATACCTTCAGCATCTGCATAAACTCAATTATGTTGGAGCTGTGTGGCATTTCACCTGCCGGG
TTCCCTTGCCAGACATGATGTCGGTTATCTGGCAAAAGCCATTTTATGTCGGCCACGAGGCAGGTGTTT
ACTGTGCATCATTCATTCACCCATGATGTTCTGGAAGTAAACCTTGGTTGTCTTCGTTGGCCAGAGG
AATTGAATGTCCGTCTGTCTGCCAGGCTGGAGTGGCACGATCTCAGCTACTGCAGCCTCCACCTCCA
GGGTTCAAGCGATTCTCCTGCCTCAGCCTCCCGGGGATCTCAAGGGCATTTCCTTTGCCACCTCACC
TTTTCATATTTGGTAACTGTATGCATTTGCCTCCAGCCCAAGATTATAAATATGAACTGATTATGATC
TGCATGTTCTCTCTTTGGGTTCAAGCATTTCTTACAGAAGAGCCACCGTGAAGTCATGGGTAAATAT
GTGTTGAATTGGTAACTCCCTCTTGGAATTTCTTGTGAATTACACAGCAATAGGGGAATCATTTAA
CTGGAGACATAATCTCAATTTGTAAAGTGTGGCCATTTTCTAACATTTTATTTTGCATACCTCCCC
TCTCTTCTCGATTGATGAACTAACAAGAGGTTAATAAAAGCCATCTCGTCATGTACAGGGACTCAT
GCCTTTAATTCCAGTGTCTTGGGAGACTAAGCAAGAACATCTTGAGGCCAGGAGTTCGACACTGTAT
AACAGATATATCTGACAGCAGTAACCTCAAGCACACCCTAAGAATGACCCTGTGGTCTAAGAAGAGTGT
GGTTCAGCGTTCCAAACGAAGAAATCCAGGAGTGGCCAACCTGGAGGTTTCTGAACCCCTGGCCTGTTC
CTTGGAATGAGGCCCTTTGTTTTGAGCTAGGCGGAGGGTGATAAGTGAAAATACTATATAAACTGCATG
CTTCTTTCAAATGATAGTGGTTTTTTTGTCCAGCCTGCCGCTCCTGGACTGCCCTGTATCTAAGTCTC
AATAAACCTATGTCTATTTTC
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_133460.3](#)

Summary: Transcriptional repressor.[UniProtKB/Swiss-Prot Function]

Locus ID: 147686

MW: 53.2