

Product datasheet for **SC214369**

Prostein (SLC45A3) (NM_033102) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Prostein (SLC45A3) (NM_033102) Human 3' UTR Clone
Symbol:	Prostein
Synonyms:	IPCA-2; IPCA-6; IPCA-8; IPCA6; PCANAP2; PCANAP6; PCANAP8; PRST
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033102
Insert Size:	1413 bp



[View online »](#)

Insert Sequence: >SC214369 3'UTR clone of NM_033102
The sequence shown below is from the reference sequence of NM_033102. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AAGAGCGACTTGGCCAAATACTCAGCGAGAAAACTTCCAGCACATTGGGGTGGAGGGCCTGCCTCACT
GGGTCCAGCTCCCCGCTCCTGTTAGCCCCATGGGGCTGCCGGGTGGCCGCCAGTTTCTGTTGCTGCC
AAAGTAATGTGGCTCTCTGCTGCCACCTGTGCTGCTGAGGTGCGTAGCTGCACAGCTGGGGCTGGGG
CGTCCCTCTCCTCTCTCCCACTCTAGGGCTGCCTGACTGGAGGCTTCCAAGGGGTTTCAGTCTG
GACTTATACAGGGAGGCCAGAAGGGCTCCATGCACTGGAATGCGGGGACTCTGCAGGTGGATTACCCAG
GCTCAGGGTTAACAGCTAGCCTCCTAGTTGAGACACACCTAGAGAAGGGTTTTTGGAGCTGAATAAAC
TCAGTCACCTGGTTTCCATCTCTAAGCCCTTAACCTGCAGCTTCGTTAATGTAGCTCTTGCATGGG
AGTTTCTAGGATGAAACACTCCTCCATGGGATTTGAACATATGAAAGTTATTTGTAGGGGAAGAGTCCT
GAGGGGCAACACACAAGAACCAGGTCCCTCAGCCACAGCACTGTCTTTTGTGATCCACCCCCCTC
TTACCTTTTATCAGGATGTGGCCTGTTGGTCTCTGTTGCCATCACAGAGACACAGGCATTTAAATAT
TTAACTTATTTATTTAACAAAGTAGAAGGGAATCCATTGCTAGCTTTTCTGTGTTGGTGTCTAATATTT
GGGTAGGGTGGGGATCCCCAACATCAGGTCCCTGAGATAGCTGGTCATTGGGCTGATCATTGCCAG
AATCTTCTCTCCTGGGGTCTGGCCCCCAAAATGCCTAACCCAGGACCTTGGAATTTCTACTCATCCC
AAATGATAATTCCAAATGCTGTTACCCAAGGTTAGGGTGTGAAGGAAGGTAGAGGTGGGGCTCAGG
TCTCAACGGCTTCCCTAACCAACCCTCTCTCTTGGCCAGCCTGGTCCCCCACTTCCACTCCCCTC
TACTCTCTAGGACTGGGCTGATGAAGGCACTGCCAAAATTTCCCCTACCCCCAACTTCCCCTACC
CCCAACTTTCCCACAGCTCCACAACCCTGTTTGGAGCTACTGCAGGACCAGAAGCACAAAGTGCGGT
TTCCCAAGCCTTTGTCCATCTCAGCCCCAGAGTATATCTGTGCTTGGGAATCTCACACAGAACTCA
GGAGCACCCCCTGCCTGAGCTAAGGGAGGTCTTATCTCTCAGGGGGGTTTAAGTGCCGTTTGCAATAA
TGTCGTCTTATTTATTTAGCGGGTGAATATTTTATACTGTAAGTGAGCAATCAGAGTATAATGTTTAT
GGTGACAAAATTAAGGCTTTCTTATATGTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_033102.3](#)

Locus ID: 85414

MW: 52