

Product datasheet for **SC214389**

MPPED2 (NM_001584) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	MPPED2 (NM_001584) Human 3' UTR Clone
Symbol:	MPPED2
Synonyms:	239FB; C11orf8
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001584
Insert Size:	1274 bp



[View online »](#)

Insert Sequence:	>SC214389 3'UTR clone of NM_001584 The sequence shown below is from the reference sequence of NM_001584. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC TTTGACCTTCCAACCCACAGGGTTCCGAAGCTCTAAATGCCCTATTGGAATGTGAGGGAAGTCTAT AAACTGCCATTTTCTAATTATAAACTTACATTCTTACTTATTACAAACCTGTGAGTCTTTTGG TAAATTGTTGGAACACAAATGATGCTAGAGGTTGTGCTTCTTATTTTATTTTAAATGGGCGATC CATTTGAAATCAGAGGAACATTGTGAATTTGTAATGACTTCTGTTTCTCAAAGGCCATGCCATTGT AAATTGTTAGTGTTCGCCAAAGGACAGCCAAGCTTTCTTTAAAAAGTGATAAAGTCTTATTTAATA TGCTTTAAGCTGAAAGAAAAAATAAGAAACAGGCAGTGTAAAAACCAACACAGATTGCACAA CTGTTTAAGAGTATTGTTGAAATATTTAATTTTCAATGTTTGTGTTGTTGTTTCTTGGAATGC TTCTTTTTCAGATGTGGTCCCAATTTATAGCAATCTTCTCAACAGAAGTAGGCATGGAAGACTTC TTTTCATACTCTCACTATAAAGAAAGCTGCATTGAGAAGAAAATGGCTGTCAATTAAGGATGGTTAA CTAGTGAGATTCTATTGTGTTATACAAGGTCTCATTGTTTGTGTTTCTTTAAATTATTTACAGCT TAAAAATACAGAAATGGAATCTGTCAAGAGCAGGTATTTATACGGTAAAAAATGAACATGCAGAC TCCTTTTCAATATGGGTTTATATATATAAGTATTTTGTGTATTATGACTACGTTAGGAGTTTATTAT TGTCAGGACAGTACAAGTCTTAAAGAGCTTGAGTTAAATCTAGGTACAGTTACAGGCATGTATAGAC TTAAATGGATGCAATGGAAGCTAACTAAATAAGGCTTAGTTGTCCTTTCTATTTAAATACCCCAAGTT GTCTTCTTACTTCTCTCCCTCTCCCATTTTGCAGTGTGTGCGATGCAATCTTCGCTAGCACAAAT ATTGTCGCTAATAGTCATTTCTGTTTTCCCATTTGTAATGCTGTTGAGCTTTATTCTATTTTATGTTAC TTTGTTAATGAAATTTAGGAAAGCAGTTGTTTCTTTAAATTTATTGTGATATTCTATATCTAGCGCCT TTATATGCAATAAAATTGCAAGATTTTAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001584.3
Summary:	This gene likely encodes a metallophosphoesterase. The encoded protein may play a role a brain development. Alternatively spliced transcript variants have been described. [provided by RefSeq, Feb 2009]
Locus ID:	744
MW:	49.3