

Product datasheet for **SC214400**

GLIPR2 (NM_022343) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GLIPR2 (NM_022343) Human 3' UTR Clone
Symbol:	GLIPR2
Synonyms:	C9orf19; GAPR-1; GAPR1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_022343
Insert Size:	1414 bp



[View online »](#)

Insert Sequence: >SC214400 3'UTR clone of NM_022343
The sequence shown below is from the reference sequence of NM_022343. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GAAGAAAACGTCTGCGCCGGAAGTACTTGTTAAATGTAATGGGAAGTGGCAGACTTAAGAACG
TGGATATGAAGTGCCTAGAACCACACAACTGGCTGTGCGTCTGTCCCTGTGGGTGTATGTGCTTGTG
TGTGTGATGCATGTGAGCGTCTCTGGCACACACTTGGACATACAGTTCTGTGTGCGCTCATTCTTAT
TACAGGAGTGAGCAAAGGAAGCATTACCCCGATGGTTACCTAGACCACGATTATTTGGATTGGGGGA
GGGGGGATCCGTTTTTTTTTTTAAATTTTGTATTTCTAAGCAAACCTCTTTGTACTTTTCTACT
TCTAATATCCATCCCTGGACTTTTTGTATTCAAATGTTTGTGATGCTGAGAAGTGAAGTTCATTTTAT
GTGATCTTCATGCGTCGTAATCTACTTTTGGTAGATAATTAAGATTATTAACCCCTCATTTAAATGTGA
CATAAAATACAGCTTTAAGCACATAAATACAAAGCAGCTCCATCAGGAACATGGAGCAGGCAGGGACT
CCATTTTACAGAATTACTGAGATTTCTCAGTTGTAAAACATGATGTCATCCTGCATGCCTCTGGAATT
CTCCAATGGGGTCGCCAAACAACAAATGGAGAAAAAAGTTTTACTTCCTTGCAATCTTCTACCTTTAA
ATAGCAAAGTACCACTACCACCACACCTCTTGCCCCCTTCCCTCTTTCTTAACTTCTGGCATTTC
GAGCTCAGCAGGTACCCCTGGTTTCTGGAGAGTTGGGCTAGGCTGAAGCTCCCCCTCCCCACCTCT
GCTAGGCAGCCAGGCCTGGTCTGGGAGACAGCCCTCACCTGCCTGGGCTCTTGCCAAAGCGGCCT
GGATGGATGAAGTCAGAGAGTGGGGTGAGGTGAGCTTACTCAGGGCCCCAGAGGAAGCCCTCAGCC
TCTGCCCTCCCCCACAGGGCGGGAGCCAGGCCTGTTCTGGCAGCTGTGGCTGCAGCTGTGCTCC
TGCTCCCTCCTGGAATGTGCGACAAGCCAAATGTTCCGGGGAGGCGGCCGGGGCAGGGGCTTAGAAG
TGCTAATATGGTTCTGTGTTTTGCCTGAAACGATACCAGGTTCCCTGAATAGCAACTTTACAAGGTCC
ATGTGGGAGGGACCAACCCAGATGCCCTGCTGAGTGCTCCTGAAACCATGGCAGCTCCATCTGTCAAGA
TGGCAGGGGCGGAGTGAGGGGGCTGCTGGCTTAACAGCAGGCATCTGGGCAGGCCAGTCTCAAAGCA
GCTCCTGAAGGTCTGTGTTGCACTGTCACCACTCTCAAGCTATGCCTCTAATTTACCAGGGATATTGA
CTAAGAAGACAATAAAATCTTTTCTTTGTGTAA
ACGCGTAAGCGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_022343.4

Locus ID: 152007

MW: 53.4