

Product datasheet for **SC214405**

MANEL (MANEAL) (NM_152496) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MANEL (MANEAL) (NM_152496) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MANEAL
ACCN: NM_152496
Insert Size: 960 bp
Insert Sequence: >SC214405 3'UTR clone of NM_152496

The sequence shown below is from the reference sequence of NM_152496. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCAAAGAGAAGGAGCAGTGGCTCATGTAGGGGGCCTGTAATGGGCGTGAGGTGCTGATGTCCTTGCC
TTGCTGGAAGATGTCACCATGTGGGTTACAGTAGGTTGTAGCCACTACTCGTCCCAGTCGCAGAGG
TCAGCAGATGGGTGTTTCTGGGTGGGCCGTCAGGCATGGGCTGTGCAACACAGAGCCCGTTCCTCAGG
CGAGTGGTGTGAGGTGCTCTGTGGTATGGGAACGGCAGAGGCTGGCAGGTGACTGAACTTAGCCAG
GGCAGTCCACATCCTGGGAACACTTTCATGTAACCCCTTCTAGTTTTGAACTCTGCAGCAGGCTGGTG
CTGTGTAGTGGCCACCCAGTCAACCCCTTGAAGGGTGTAGGGTCTGGGCTCGCATGCACCCTGCTC
CCTTCTGCCAGCTGCATCCTCCTCAGGCCTCCTTCCACATCATCTGCTTCTCTAAGTTAGGGAA
CCATATTTGAGACTTTCAAAAAGGGAACCTTCTAGGTTTAAAGCTCCTCCAGTAGATTCTGAACCAAT
TATCAGGAAATCATCCTGAGCACTCACAGGTTCATTTAACACTCACTCATCAAGCACCTTCTATGTGC
TAGGTGCTGGGGAAAACCTTGACCAAGGAAGAGTTCCTGTCTGAAGCTTCCAGGACCCAGCTTCCCTT
TCTGGTGTGGCCTTGTAGCTAGTGCCTGGGCACAGGTGTTTTTCTTTTGCAGTTTTACCTAGTGCTGG
GAGTTCAGTTCCTTTTCTCTAGAAAAATACCTCTGTGCTCCAGAGCCTAATTTTTCCAGATGCATAT
TTAGCTCTAGGGAGAGGACTAGGAGGAAATCCCCCTCCCTTTAGCTGCCTGAACTGACTGAGGCCACT
CACTAGAGCCATGTTCAGTGCTACTGTGATTAGTAGTAATTAATATGAACTGGTATTCTCAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_152496.3](#)

Locus ID: 149175

MW: 35.5