

Product datasheet for **SC214407**

MANEL (MANEAL) (NM_001113482) Human 3' UTR Clone

Product data:

Product Type: 3' UTR Clones
Product Name: MANEL (MANEAL) (NM_001113482) Human 3' UTR Clone
Vector: pMirTarget (PS100062)
Symbol: MANEAL
ACCN: NM_001113482
Insert Size: 960 bp
Insert Sequence: >SC214407 3'UTR clone of NM_001113482

The sequence shown below is from the reference sequence of NM_001113482. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCCGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAACGATCGCC
ATCAAAGAGAAGGAGCAGTGGCTCATGTAGGGGCCTGTAATGGGCGTGAGGTGCTGATGTCCTTGCC
TTGCTGGAAGATGTCACCATGTGGGTTACAGTAGGTTGTAGCCACTCACTCGTCCCAGTCGCAGAGG
TCAGCAGATGGGTGTTTCTGGGTGGGCCGTCAGGCATGGGCCGTGCAACACAGAGCCCGTTCCTCAGG
CGAGTGGTGTGAGGTGCTCTGTGGTATGGGAACGGCAGAGGCTGGCAGGTGACTGAACTTAGCCAG
GGCAGTCCACATCCTGGGAACACTTTCATGTAACCCCTTCTAGTTTTGAACTCTGCAGCAGGCTGGTG
CTGTGTAGTGGCCACCCAGTCAACCCCTTGAAGGGTGTAGGGTCTGGGCTCGCATGCACCCCTGCTC
CCTTCTGCCAGCTGCATCCTCCTCAGGCCTCCTTCCACATCATCTGCTTCTCTAAGTTAGGGAA
CCATATTTGAGACTTTCAAAAAGGGAACCTTCTAGGTTTAAAGCTCCTCCAGTAGATTCTGAACCCAAT
TATCAGGAAATCATCCTGAGCACTCACAGGTTTCAATTAACACTCACTCATCAAGCACCTTCTATGTGC
TAGGTGCTGGGAAAACCTTGACCAAGGAAGAGTTCCTGTCTGAAGCTTCCAGGACCCAGCTTCCCTTT
TCTGGTGTGGCCTTGTAGCTAGTGCCTGGGCACAGGTGTTTTTCTTTTGCAGTTTTACCTAGTGCTGG
GAGTTCAGTTCCTTTTCTCTAGAAAAATACCTCTGTGCTCCAGAGCCTAATTTTTCCAGATGCATAT
TTAGCTCTAGGGAGAGGACTAGGAGGAAAATCCCCCTCCCTTTAGCTGCCTGAACTGACTGAGGCCACT
CACTAGAGCCATGTTCAAGTACTGTGATTAGTAGTAATTAATATGAACTGGTATTCTCAA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).



[View online »](#)

Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001113482.2</u>
Locus ID:	149175
MW:	35.5