

Product datasheet for **SC214427**

GCC1 (NM_024523) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	GCC1 (NM_024523) Human 3' UTR Clone
Symbol:	GCC1
Synonyms:	GCC1P; GCC88
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_024523
Insert Size:	1416 bp



[View online »](#)

Insert Sequence: >SC214427 3'UTR clone of NM_024523
The sequence shown below is from the reference sequence of NM_024523. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GCCAGCTGGTGGCCTTCTGGCAAGAGATGATGCCATTTTCATGAATTCCTGGAATTTCTGTGCCTCTTT
TATGTGGGAAAGGGCAGGCTCTGGTTTCTGCTGCCTCTTCTTACATTGTCATTATTTCTTCACTGT
GTTGAGCTGGGAGGAGTCTTCAATGTGGGAAGGGATTTTCTGCCAAATGCCATTAATGCTGCTTTATCC
TTAGGCCCTAGAGATACTGAAAGCTGGGGCAGTATCTTCTGATGGTGTGCATGTGAGGAGAAGATGA
GGTTAGGACTGAGAAGTGCAGAAGTTGGAACAGTGGTAAGGCTGTTTTAAATAAGATGTTTTGTTTA
ATAATATGCTCCTGGCACAAGCTAGGAGTAAATGTGACTCCAAGGGAGTTCAGTTAATCTCTGAAAT
GCACAAAACCTAGCTATTTTCTCCCTCTCATCAGTCTGAGTCTGGTCCATTGCTACCCCAATTCTCT
GGGGACATAAAACCAGGCTGGAAAGGGACCAGGAAGTTTGAAATAGTGACATATCATCCACTAGTCCCA
AGGGCTAAGGAATAGTGAGTTTATTCTGGAAGGAAGCTGGGAAGCTTAGTCTAATTAGTGCCTGGGGATG
ACCTATGCAATCACACCGCTTATGACCATCCTAGAGAGGGCCCTGAGCACCAGCTTGATCTTAGGGATT
TCCAAAGTAACCTGCTTTTTGCCTGGATAGGGTTAAAATAGACCTTTCTTGCCTATCCTTGCCTTAACC
TATCTGCCTGAGGTTGGCCTGAGATTGTGAGTCAACGACTTTGCTATCTTTCTCAGTGTTGAACCTTT
CATTAAGAAATAAAGTCCTAGCTTCTTACAGAGAGGGGTCCAAATGGTGAATGCTCATCTGCCTGGAT
TCAAGGAATTAGCTCAGAGATTGGCCCCAGCTTTTCTGCCTTTGTAGGGACAGCAAAAGGGGAAATTT
TGCTGCAGAAAATTCCAAAGATTGCTGTAGCTCTCACAGGGAAGTGGTAAAGATCAGCTAAACCTGGG
TTGGGGTGCTTTCTGCCAGTGGGTCTTGGCATAAGTAGATTAATCCTGCTCTTTAAGAAAAGGCAAC
TTATTCAGGCAGTCTGGAAAGGGGTCTCAGAAAACCTAGTTTCTTTATTCTTCTTTTCTCCCACT
ACTGTTACTGGTTATAGAGGTCTTTGGACTCTAAAGACCAATGTTTGGCCACTAACTGGACTAATATGT
ATCTTTCTGTGATTTTCATCATAGAGGTCTGTTTTGTGAGGGTTTGGGGTGCAGAAAACCTTGATTAAT
CTTAATGGGAGGCTGGGTGACCTGGATTATCTACAGTGAGCAGACTTAAATGGAACAGAAGTTTATGTG
TCCAAATGATGGAATCATTAAACCTGAGTGACTTGA
ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_024523.6](#)

Summary: The protein encoded by this gene is a peripheral membrane protein. It is sensitive to brefeldin A. This encoded protein contains a GRIP domain which is thought to be used in targeting. It may play a role in the organization of trans-Golgi network subcompartment involved with membrane transport. [provided by RefSeq, Jul 2008]

Locus ID: 79571

MW: 52.9