

Product datasheet for **SC214439**

ERI2 (NM_001142725) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ERI2
Synonyms:	EXOD1; ZGRF5
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142725
Insert Size:	1414 bp



Insert Sequence:

>SC214439 3'UTR clone of NM_001142725

The sequence shown below is from the reference sequence of NM_001142725. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TTGAGACTCAGGCTTCAATGAGGAATGTATAACCTTTTCATGTATGAATCCTAATTGTTCTTGAATTT
CCAAACATGAGTATTCTGATAACATCTTACACTATTTTATTTTATTTTATATTAATGTTTACTGCA
TTTACCATACTTTCTTCCAGGAACCAAGTGAATCTTTTGTAAACAGTATTTATTAGATTACATATAA
TAAGCAACAATAGTATATGATCCTTCTTTTATAGTTAGACCTTGAAACAAGTCTAAACTATTTTAAAT
TATTTGCAAAAAATGTGACTAATAAAAAATCCCTAGTCTATTTTCAGGTAGATTCTTGTACTTTGAGATT
ATCATAATCTTCTAACTGATATTTTATTTTCTAAAACTTTTATAAGTTCTTTAGTTTTATTTCAG
GCTTGAACATTGGTTCTGATTGTATATCAACATATCAGAAGAATTAGTCAACTAAAGTTTGAAACCT
ATGTAGCAGAGATCTGGTTCAGAACTGCTTCTATACTATACATTTAAAAATTAATCAGCAATTGATGG
TATTATTTTATATTAAGAAATGATAGGAACATTGACACTAATGTGGAGGAATTTGCCAAATATTTTAT
TTCAAAATTTAGCCTCTGATATCATAACCACAAGATATTTTAAATCAACAGTTTCTAAGCATACTACCA
CAGTTTTTATGACCTCTCCATCATCTTGTGGTTTCTAGAGATTGACAGATACTATGTTTGTTTTATAGA
TATGGTAATTAATGAAACAACTTTAAATTTGCTTCCATTCTTCTTCTCAGTTCATTTCTTTTGTCT
TCCCAGTGATAGTCTTTGGCAGCTCTTGAATAAATTCTACCTGTTGGCATCAAGGGGAAGAAATTAGC
CTGGAATTTTATCTTTACATAAAGCCAATTATATAGCAGTTTGATTATGAATTCTAGGTGGAAGAAT
TAGAATAGTCTTTGGAATGGGTTTGGGTAAACGCCAAAATAGGGAAAAGGGTATTGTTAGGAAAGTTG
TCTTGATATTAAGTATTTACTGTGAGTCCAGTAAGCCAAAGGACCAACAAAATAGGTTCAGGGTTTCA
TTTTAATCTCAAAGATTCTGATTCTTACCCATATGTCATCAGTTCAGTTTGGGACATGTGGTGCAGCAT
GTGTGAAGAAAATATTCATTAATAATGTAAAAATAAACACTTGAAATTGTCCACAACACTGAGCAAATAT
TCGTTATAATTAGGATGCCTACCTTCTGGGATATTTGTAAGGTGCTGTAGTTTTTTAACATGCTCCT
GAATCTCCTTTATAGTTGTTCTTGATCATGTGACTTGTAAATCAGGATTTAGAACGACAAAAGCCTTTA
CTACCTAAAAATAAATATTTGAAACACCAGGAAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
    
```

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_001142725.2

Locus ID:

112479

MW:

55.1