

Product datasheet for **SC214445**

EVI1 (MECOM) (NM_001105078) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	EVI1 (MECOM) (NM_001105078) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	MECOM
Synonyms:	AML1-EVI-1; EVI1; KMT8E; MDS1; MDS1-EVI1; PRDM3; RUSAT2
ACCN:	NM_001105078
Insert Size:	1432 bp



[View online »](#)

Insert Sequence: >SC214445 3'UTR clone of NM_001105078
 The sequence shown below is from the reference sequence of NM_001105078. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGTGCTATCCAGTCCATAAGCCACGTATGACGTTATCAAGTTGACCAGAGTGGGACCAAGTCCAACAG
TAGCATGGCTCTTTCATATAGGACTATTTACAAGACTGCTGAGCAGAATGCCTTATAAACCTGCAGGGT
CACTCATCTAAAGTCTAGTGACCTAAACTGAATGATTTAAAAAGAAAAGAAAAGAAAAGAAAAGAACTAT
TTATTCTCGATATTTGTTTTGCACAGCAAAGGCAGCTGCTGACTTCTGGAAGATCAATCAATGCGACT
TAAAGTGATTCAGTGAACAACAAAACTTGGTGGGCTGAAGGCATCTCCAGTTTACCCACCTTAGGG
TATGGGTGGGTGAGAAGGCAGTTGAGATGGCAGCATTGATATGAATGAACACTCCATAGAACTGAAT
TCTCTTTTGTACAAGATCACCTGACATGATTGGGAACAGTTGCTTTTAATTACAGATTTAATTTTTTTC
TTCGTTAAAGTTTTATGTAATTTAACCCCTTGAAGACAGAAGTAGTTGGATGAAATGCACAGTCAATTA
TTATAGAAACTGATAACAGGGAGTACTTGTCCCCCTTTTGCCTTCTTAAGTACATTGTTTAAAAGTAACTAG
GGAAAAAGGGTATGTGTATATTGTAACACTATGGATGTTAACACTCAAAGAGGTTAAGTCAGTGAAGTAA
CCTATTATCACCAGTACCGCTGTACCACTAATAAATTGTTTGCCAAATCCTTGTAAACATCTTAAT
TTAGACAATCATGTCAGTGTGTTTTAATGTTTATTTTTTGTGTGTGTTGCGTGTATCATGTATTTATT
TGTTGGCAAACTATTGTTTGTGATTAATAAGCACTGTTCCAGTCAGCCACTACTTTATGACGCTGA
GGCACACCCCTTCCGAATTTCAAGGACCAAGGTGACCCGACCTGTGTATGAGAGTGCCAAATGGTGTT
TGGCTTTTCTTAACATTCCTTTTGTGTTGTTTGTGTTTGTGTTTCTTCTTAATGAACATAAATACGAATAG
ATGCAACTTAGTTTTTGTAACTGAAATCGAATCAATTGTATAAACGATTATAATTTCTTTCATGGAA
GCATGATTCCTTCTGATTAATAAACTGTAACCATATTTTATGCTGTTGTCTGCAAGCTTGTGCGATGTT
ATGTTTCATGTTAATCCTATTTGTAATAAAGTGTCCCAACCTTATGTTAAAAGAGAGAAGTAAATAA
CAGACTGTATTAGTTATTTGCCCTTTATTGAGGAACAGATTGTTTTCTTTTGTGTTAATCTCA
TTTTGAAATAATCAGCAAGTTGAGGTAATTTCTCAAATGCTTTGTACAATATAAACTGTTATGCCCTT
CAGTGCATTACTATGGGAGGAGCAACTAAAAATAAAGACTTACAAAAAGGA
ACGCGTAAAGCGCCGCGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001105078.4](#)

Summary: The protein encoded by this gene is a transcriptional regulator and oncoprotein that may be involved in hematopoiesis, apoptosis, development, and cell differentiation and proliferation. The encoded protein can interact with CTBP1, SMAD3, CREBBP, KAT2B, MAPK8, and MAPK9. This gene can undergo translocation with the AML1 gene, resulting in overexpression of this gene and the onset of leukemia. Several transcript variants encoding a few different isoforms have been found for this gene. [provided by RefSeq, Mar 2011]

Locus ID: 2122

MW: 55.5