

Product datasheet for **SC214461**

STAG3L4 (NM_022906) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	STAG3L4
Synonyms:	FLJ13195
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_022906
Insert Size:	1445 bp



Insert Sequence:

>SC214461 3'UTR clone of NM_022906

The sequence shown below is from the reference sequence of NM_022906. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CTGTTGGAGAAACACAAAGAGCTGCATTAACTACTACATTGTCAAAGACAAAGACGAACATGAGACATGC
GATATGAGTTTTAAAGGAGAGAGAGCAAAAAACAATACAGCTATCCAAGGAGCAAGTGCTGAGCTT
GTGAAAAACAGCTTCTATTTTATTAATAATAGTGTTCTCGCAGGGATGTCTGAAGTACTCTGCAAGG
CTAATTAGGACAGATCAATGTACACTTACCTCTCCAGAAGGCTGAAAGGTGTGTAGTGACTATGTGTCA
TTAAAAATTTACTTGACCATGGGCAACAGATGTCCATGAAAGAATGACCTGGTGAGTGCTCTAAGTTT
TCATTTTACTCATGTATCCTCAATGCATCATAAATGCCATAGACACCATTACTATCACGTAGGCTCAT
AGCAGCTGATGTGCTGTATATGCTTGAGAGGGTGGTATCAATTAATATATTATCTATAAAGCCATACT
TGCAATTACTACCACTTCACAATTTATCATTAAACGCATTTGGCAGATTTTTTTCTCCAAAGTAACAT
TTATTTAAAGGTTTATCAAAATCTTAGAACAGCATTATTATAATGTGATTTTGTAAATTGACTAAA
TCAGTGGCTCTTAAAAACACCTCATCTCAGTCCCGCTCAGTATACACCACAGTGCCATTTGGGAAAA
GAAATACATATTAGAAATGTTTCATTCAGACCTCGTCATGGAAGATCATATTAATGTAATCAGGTGG
TACGTTTACATGCTCATTAGTAATCCATCATTAGTCCAATGCTCACTATTAGGAATAAATGTACAACG
AGTTATTCAGAATACCCAGAAAGAGTATTTCTTATTTAAACATTCAGGACAAGCTGACTGAAATAAGTG
CACACCTTGAATTTGAATTTCAAGGCACATCACTAAAGTCAATCAGTTATGTCAGGACACTCTTCTAT
AATAATCTCAATTATAAAATCTCGCTCTCCCTAACACAGGAGTTGGAGGGAGTTCCAATATCTGAAT
CTTTCTAATTCAGTGAATGCAAAAAATCTTGATGAAGTGTGATTTAAATTTTCATGCCTTAATACA
TAGCATGCCACTAAAACACATTAGTTATGCTTTATGCTGCTTAATCATACGTCAAAATCTGTAAATC
TGCTCTTGATGCTAAGGTAATTTAACTGCAATTAATCAATTAAGAAAAATAAAATGTAATAAATGCCA
AATAATGTAATCTTTTAAACAGACCACTATCTACACATTTATCCCCAAATCTTATTAGCTGACAC
CATATAAGGAAAAACATAATAAAAGTAGACTTCACTGAGCCTTTCTCCATTAGTAGCTTGCTCTAAAC
TTACAAATGTTTTAGAGAATTCTTTGGCATGATCTGACTTTCCAATAAAATGTGAGCAACTTGA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_022906.2

Locus ID:

64940

MW:

56.4