

Product datasheet for **SC214465**

ZNF781 (NM_152605) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF781 (NM_152605) Human 3' UTR Clone
Symbol:	ZNF781
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_152605
Insert Size:	1444 bp



[View online »](#)

Insert Sequence:	>SC214465 3'UTR clone of NM_152605 The sequence shown below is from the reference sequence of NM_152605. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC AGAAGGCCTTCAAGCATAATGCACATCTAATGAGACATCAGAGAATTCACAGTGGGAGAAACCTATG AATGTAATGAATGTGGGAATGCCTTCATTTGTTTTTTGATTTTAAATAGACATCAGAGAATTCACAAAG GAGAGAAACCTATCAGTGTAAAGAATGTGGGAAGTTTTTAATAATTATCCAACCTTACTCAACATC AAAGAATTCATACCAGTGTGAAACCTATGAATGTAAGGAATGTGTAAGACCTTTACTTGGCTTTCAA ACTTGAATAAACATCAGACAATTCACACTGGTGACAAACCTATGAATGTAAGATTGTAAAAAGCTT TTCTTTGTTATTCAACATTTATTCAACACTATAGAACTCATACTAATGAGAAACCTATGAATGTCAGGA ATGCAGCAAGGCCTTTAAGCAGAGAGCCCATCTAGTCAACATCAAAATATTCATACTGGTGAGAAACC CTATGAGTTTAAGGAATGTGGGAAGGCCTTTGCTTGTTCCTACCTTAATAGAAATCAGAGAATTCACA CTAGTGAGAAACCGTATCAGTGTAAAAATGTAGAAAATCCTTTATTGATTGTTCAACTCTTAACACCA GAGAATTCACACTGGTGAGAAACCTTTGAATGTCAGGAATGTAGGAAGGCCTTTAGGCAGAGGGCGCA TCTTACTCAACATCAAAGAATTCCTACTGGTGAGAAACCTATGAATGTAATGAATGTGAGAAGGCCTC TTCTTGTCTTTTATTTAACAATCCTTATCTAACATCACCAAATTCACATTGCAGAGAACTTCAGAAA TATAAGAAATGTGGGAACACTTTTACCTGTCCCTCAGCTATACTTTTCATAACAGAATTTGTACAGG TGAGAAGCCTAATGAAAGCCGATATGTGTAAGATCTACTCTCACAATATACTCAACATTATATAATG TATAGTGTGATTAATATAGGAAATCCTTGATTTTGAATTATCAATTATGGAATAGCAACAATTTGTA AGTGAAGAAAATCCAGTGACTGTAAATTTTTTACCTCATAATTCACAATTTATTAAGCATCAGAATTT ATATGGAATAATAGAATATCAATTATTAATAAAAAACATTACCACTAAAAATTTTATTGACTTATAGAG AATAGCATATAGAGAATAATCATGTTGACTTTGTGAATGTAGAAAAGCTTTTATCTATTACAGAACTCT TTTTCAGCATCATTATACTTTTACCAGTTTCCACGCTGTGTTTATTACCAAGTAATTTAACTGTATAC TTCAGATTACAAATTTCTTACATTGTGATAAGGTTATTGTGAAGATTAATGAGTTTATCCAAA ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_152605.4
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	163115
MW:	56.8